

Maxenon™ Xi300



NEVER COMPROMISE  
START WITH YOUR VISION

Extreme and dependable  
xenon illumination exceeding  
the standard



We perfected  
the 300-watt xenon  
light source

Since xenon's advent in surgery during the late 1980s many 300-watt xenon light sources presented problems during time-consuming cases. We perfected the 300-watt xenon light source.

Our Maxenon™ Xi300 is so bright many surgeons turn it down – confidently knowing they have higher gears to shift into when the going gets tough. More important are the benchmark features such as a true Hot Start and the technologically advanced fan we build into our 300-watt xenon system that prevent failed start ups and overheating. The result is dependable, extreme white xenon light through the longest and most difficult cases. Our signature five-port rotating turret enables use with ACMI, Wolf, Storz, Olympus and BFW cables making Maxenon™ Xi300 the most versatile xenon light source in its class.



NEVER COMPROMISE  
START WITH YOUR VISION



The FOS 4X Xenon System combines our advanced High Beam TriLens fiber optic headlight with an intensely bright, dependable and durable 300-watt xenon light.

## Specifications

|                      |                                                                                                                                                                                                                                    |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Lamp</b>          | 300-Watt Xenon<br>Color Temp. 6,000 <sup>o</sup> Kelvin<br>650-Hour Lamp Life                                                                                                                                                      |
| <b>Dimensions</b>    | 12.6" L x 11.5" W x 5.25" H<br>(30.5mm x 29mm x 13.5mm)                                                                                                                                                                            |
| <b>Weight</b>        | 11 lbs. (5 kg)                                                                                                                                                                                                                     |
| <b>Power</b>         | 100-240V<br>50/60 Hz                                                                                                                                                                                                               |
| <b>Compliance</b>    | ETL Certified to UL 2601<br>CETL Certified to CAN/CSA<br>C22.2 No. 601.1<br>CE Certified                                                                                                                                           |
| <b>Product Codes</b> | MID3010<br>Maxenon™ Xi300 Light Source<br><br>FOS 4x<br>Maxenon™ Xi300 Light Source<br>With High Beam TriLens Headlight<br>7' cable and Mobile Stand<br><br>MID3020<br>Replacement Lamp and Housing<br><br>MID2004<br>Mobile Stand |



1800 Williamson Court  
Louisville, Kentucky 40223 USA  
www.dremed.com  
**(502) 244-4444**  
**(800) 477-2006**