

Maxenon™ PowerPlus



NEVER COMPROMISE
START WITH YOUR VISION

Remarkable quality and powerful
xenon light within your means



PowerPlus technology
affords a degree of
brilliance you will be
amazed to see

Well suited for hospital operating rooms or outpatient surgery centers, the Maxenon™ PowerPlus will drastically improve headlight intensity beyond expectation. Our innovative PowerPlus technology affords a degree of brilliance you will be amazed to see in a 180-watt xenon light source. Our signature five-port rotating turret enables use with ACMI, Wolf, Storz, Olympus and BFW cables making Maxenon™ PowerPlus the most versatile xenon light source in its class.



NEVER COMPROMISE
START WITH YOUR VISION



The FOS 3 Xenon System combines our advanced High Beam TriLens fiber optic headlight with the Maxenon™ PowerPlus, resulting in a 180-watt xenon system that is reasonably priced and beyond compare.

Specifications

Lamp	180-Watt Xenon Color Temp. 6,000° Kelvin 650-Hour Lamp Life
Dimensions	12" x 6 3/4" x 6" (205 mm x 170mmx152 mm)
Weight	7.75 lbs. (3.5 kg.)
Power	100-240V - 50/60 Hz
Compliance	ETL Certified to UL 2601 CETL Certified to CAN/USA C22.2 No. 601.1 CE Certified
Product Codes	MID 1810 Maxenon PowerPlus FOS 3 Maxenon PowerPlus Light Source with High Beam TriLens Headlight, 7' Cable and Mobile Stand MID 1805 Replacement Lamp and Housing Drawer MID 2004 Mobile Stand MID 4005 Carrying Case for Light Source and Headlight



BFW, INC
2307 RIVER ROAD, STE. 103
LOUISVILLE, KY USA 40206
WWW.BFWINC.COM
502.899.1808
800.717.4673