


Patient monitoring solutions for your medical facility



Transport : Emergency



Bedside : Vital Signs



O.R. : Critical Care



Central Station Monitoring



When monitoring
multiple parameters
 & **multiple patients,**
 we have you covered.

DRE Envoy TS
 Monitor patients throughout your facility from a single location

The configurable DRE Envoy TS Central Monitoring Station includes wireless ECG transmitters. No other system provides the high level of flexibility offered by the DRE Envoy TS. Each channel is independent, allowing the use of stationary, wireless transport and ECG telemetry transmitters on a single system.

Parameters: ECG monitoring using portable transmitters, or get complete integration with DRE-brand monitors for ECG, respiration, temperature, SpO₂ and NIBP.



DRE Waveline EZ

Portable Patient Monitor with an intuitive touch-screen

At only 6 pounds, the Waveline EZ is one of the most mobile multi-parameter monitors on the market. It features an 8" high-resolution touch-screen that provides immediate operation.

Parameters: ECG, resp., SpO₂, NIBP and temp. Available with EtCO₂ and a printer.

Designed for: Triage, emergency environments, ambulance use, and patient transport.

DRE Waveline Touch

Lightweight Patient Monitor with a large touch-screen

The fully-functional Waveline Touch features a 12.1" color TFT display with a touch-screen interface. It features multi-lead ECG monitoring, ST segment monitoring and 72 hours trend data analysis.

Parameters: ECG, resp., SpO₂, NIBP, heart rate and temp. Available with dual IBP, EtCO₂, and a built-in recorder.

Designed for: Pre-operation and post-operation, ER, transport, CCU, and conscious sedation rooms.

DRE Waveline Pro

O.R. Monitor available with EtCO₂ and anesthesia gas monitoring

The lightweight Waveline Pro displays eight colorful waveforms on a crisp 15" high-resolution screen. It works with EMR systems, and the optional agent bench is part of main unit—making it a very compact anesthesia monitor.

Parameters: ECG, resp., dual temp., SpO₂ and NIBP. Available with EtCO₂, IBP, anesthetic gas monitoring, and a printer.

Designed for: General surgery, ICU, and anesthetic agent monitoring.

Envoy Telemetry Transmitter

Portable ECG transmitter for Envoy TS Telemetry System

Flexible, durable DRE Envoy TS telemetry transmitters feature user-selectable lead wire configurations and provide an extended battery life.

1-800-477-2006 • www.dremed.com





Monitors | Anesthesia | Electrosurgical | Surgery Lights | Surgery Tables and much more



DRE Waveline EZ
Portable Patient Monitor



DRE Waveline Touch
Patient Monitor



DRE Waveline Pro
O.R. Monitor



DRE Envoy TS
Central Monitoring Station

1-800-477-2006 • www.dremed.com

