


Thru-the-Lens Headlight Video Imaging System



NEVER COMPROMISE
START WITH YOUR VISION



The Thru-the-Lens Headlight Video Imaging System is available with its own tower cart including a medical grade flat screen monitor and DVD recorder.

Additionally, our camera is custom engineered for use with a surgical headlight – not a refitted endoscopic or pre-built surveillance camera. This is how we offer incomparable single-chip image quality. With DV (IEEE 1394), S-VHS and BNC capability, our system offers an extended depth of focus. The image remains in focus regardless of inward or outward movement from the field of surgery.

The Thru-the-Lens Headlight Video Camera System is available in 180 or 300-Watt Xenon Formats.

Specifications

The Thru-the-Lens Headlight Video Camera system includes:

- Thru-the-Lens Headlight with 10' Fiber Optic Cable
- Camera Control Unit
- 180-Watt PowerPlus Xenon or Maxenon™ Xi 300-Watt Light Source
- All necessary Video Cables and Power Cords
- Mobile Stand

Product Codes

- MID 7910
Thru-the-Lens Headlight Video Camera System (NTSC)
With Maxenon™ PowerPlus 180-Watt Xenon Light Source
- MID 7915
Thru-the-Lens Headlight Video Camera System (PAL)
With Maxenon™ PowerPlus 180-Watt Xenon Light Source
- MID 7930
Thru-the-Lens Headlight Video Camera System (NTSC)
With Maxenon™ Xi 300-Watt Xenon Light Source
- MID 7935
Thru-the-Lens Headlight Video Camera System (PAL)
With Maxenon™ Xi 300-Watt Xenon Light Source
- MID 8610
10' Replacement Fiber Optic Cable
- BFW 1505
Medical Grade DVD Recorder
- BFW 1516
Flat Screen Monitor
- BFW 1520
Tower Cart



BFW, INC
2307 RIVER ROAD, STE. 103
LOUISVILLE, KY USA 40206
WWW.BFWINC.COM
502.899.1808
800.717.4673