




Available from  **DRE**

Westar Medical Products, Inc.
OS III Oral Surgery Chair



ELEGANCE & DISTINCTION Westar introduces affordable luxury. Always famous for reliability and efficiency, Westar adds beauty and comfort to the dental operator. The latest technology is combined with the skill of craftsmen to create hallmark products which move you to the forefront of modern operatory design.



Toll-free **1-800-477-2006** Local **502-244-4444** Online www.dremed.com
DRE, Inc. 1800 Williamson Court Louisville, KY 40223 USA

OSIII Oral Surgery Patient Chair

- ❖ OSIII with Asepsis upholstery
2000-080-2

Order#

Standard Features:

Electro-mechanical cantilever system
 Foot control with two programmable preset positions
 and automatic exit (return)
 Wheelchair accessible double articulating headrest
 with lock-down for surgery feature
 Left & Right IV arms with one interchangeable pillow
 HiGlos powdercoat inside and out
 Thermal overload protection
 Synchronized seat and back tilt
 Low voltage
 Three year warranty
 Height to seat: 13" to 31"

Factory Options:

Eight way articulating I. V. arm
 Touchpad with 5 programmable preset
 positions and automatic exit
 Headrest options
 Additional I. V. pillow only
 3" Height adjustment spacer...2000-119
 Leg/Body Strap set.....2000-223
 Mobile Base.....7000-371
 Ultraleather upcharge (asepsis only)1005



Shown: OSIII with optional eight way I. V. arm and Chinner Headrest

2001 Consultation Chair

- 2001 with Asepsis upholstery ❖
- 2001 with Sling upholstery & Trim ❖
- 2001 with Plush upholstery* & Trim ❖

Order# 2000-401

*Plush: 2000-401S

2000-401T

Standard Features:

Electro-mechanical cantilever system
 Foot control with two programmable preset positions
 and automatic exit (return)
 Wheelchair accessible double articulating headrest
 Left/Right patient access with arms that swing up
 HiGlos powdercoat inside and out
 Thermal overload protection
 60 degree swivel and brake
 Synchronized seat and back tilt
 7 way foot control
 Low voltage
 Three year warranty
 Quick change upholstery

Factory Options:

Touchpad with 5 programmable preset positions
 and automatic exit (return)
 Non-articulating headrest in lieu of articulating
 Covered arms for Asepsis upholstery (set)..7000-096-1
 (*Plush upholstery done in Ultraleather only)



Shown: 2001 Electro-mechanical patient chair with asepsis upholstery and optional covered arms

Factory Options for Patient Chairs



- ❖ Non-Articulating Headrest
in lieu of double articulating
7000-028
- ❖ Ponytail headrest upgrade
7000-609



- ❖ Chinner headrest
7000-608



- ❖ Eight way articulating I. V. arm
2000-220



- ❖ Touchpad with 5 programmable
preset positions and automatic exit