



Certified Refurbished

The durable, flexible choice for your operating room

The Amsco 3080 features a radiolucent top and can be adjusted to Trendelenburg, back, side, foot and flex positions. It comes standard with a remote control for easy patient positioning.

- **Amsco 3080 RC** can support up to 400 lbs.
- **Amsco 3080 RL** offers return to level option, and can support up to 400 lbs.
- **Amsco 3080 SP** offers return to level option, and can support up to 500 lbs.

► A Quick Look

- Simple patient positioning
- Supports up to 500 lbs.
- Radiolucent top
- Durable, flexible and reliable
- Fluoroscopic capabilities
- Returns to level with a single command

Product Specifications

Load limit	up to 500 lbs
Maximum height	44"
Minimum height	27"
Trendelenburg	25°
Reverse trendelenburg	25°
Lateral tilt	18°