

## High frequency desiccator with 60 watts of cutting power and coag capabilities



It's half high frequency desiccator, half Bovie Generator. It's the Aaron 950, and it can be utilized for almost any cutting procedure in the doctor's office, including LLETZ/LEEP.

### **Three-Button Pencil**

A masterpiece of engineering, the Aaron 950 features a three-button pencil that allows you to adjust the power output up and down. You can also use the pencil to activate the unit, even when you are in foot control mode.

### **Designed For Years of Trouble-Free Use**

The Aaron 950 is designed and engineered to be the most reliable and durable generator available today. The unit's life is prolonged by incorporating the latest in chip technology and reducing the number of internal parts, thereby generating less internal heat. Cleanup is simplified by its non-absorbing, nonstaining, high-impact ABS enclosure.

#### **► A Quick Look**

- High frequency desiccator with cut
- Allows 60 watts of cutting power and coag capabilities
- Includes a three-button pencil that allows you to adjust the power output up and down, and activate the unit
- Engineered for reliability and durability
- Features the latest in chip technology
- Easy to clean with a non-absorbing, nonstaining, high-impact ABS enclosure.

### Technical Specifications

#### The Generator

The left half of the unit is a Bovie generator. It is a floating ground unit that requires a return electrode pad to be used. This isolated power feature helps to ensure your safety as well as that of the patient.

Simplicity is achieved using easy-touch, color-coded membrane pads, and by placing all user controls and outputs on the unit's compact front panel. The ergonomic arrangement of commonly used operating functions, plus a large, comfortable rotary dial, provide quick operating-mode selections and smooth power adjusting with the option of adjusting the power from the handpiece.

#### Digital Error Detection

Digital error detection means unsurpassed safety for you and your patient. The Aaron 950 continuously monitors every aspect of the unit's output. At the sign of any problem, the machine instantly disables the output and displays the appropriate error code on the display.

#### Six Presets

The six available presets labeled "A-F" allow the user to record often used settings for one touch recall. Each preset stores the mode selection and power setting for one operating mode.

#### Three Modes of Operation

The Bovie Generator offers you up to 60 watts in Cut mode, Blended Cut mode, and Pinpoint Coagulation mode.

#### Dimensions & Weight

Width ..... 26 cm (10.25")  
Height ..... 17.8 cm (7")  
Depth ..... 11.4 cm (4.5")  
Weight ..... < 3.7 kg (< 8 lbs.)

#### Output Power

##### Generator

Cut ..... 60 watts  
Blend ..... 60 watts  
Coagulation ..... 60 watts

##### High Frequency Desiccator

Fulguration ..... 35 watts  
Bipolar ..... 30 watts

##### Line Voltage

A950 ..... 120 VAC  $\pm$  10%  
A950-220 ..... 230-240 VAC

Output Frequency ..... 350-800 kHz

Line Frequency ..... 50-60 Hz

#### Available Accessories Include:

- Handpiece (A801-EU)
- Handpiece Drapes (A910)
- Footswitch (A1203)
- Footswitch, Waterproof (A1203W)
- A1205A Adapter for connecting Footswitch Pencil (ESP7)
- Bipolar Cord (A827-EU)
- A1202C Connecting Cord for Dispersive Electrode (Grounding Pad)
- A1202 Solid Dispersive Electrode (Grounding Pad)
- A812-C Mobile Stand with Bottom Smoke Evacuator Tray, Top Instrument Tray and Clamp, and 12" Height Extender
- A812 Mobile Stand
- A812-BT Bottom Tray for A812 Mobile Stand
- A808-T Top Instrument Tray and Clamp for A812 Mobile Stand
- A8121-E 12" Height Extender for A812 Mobile Stand
- A813 Table Stand
- A837 Wall Mount Kit