

FAMILY PRACTICE TABLE

For comfort, economy, overall functional and storage, there's nothing better than Clinton's Family Practice Table. This versatile table has storage access from either side. Plus, its lower treatment table height could eliminate the need for a step stool, and makes it ideal for family practice, clinic or wherever a cost effective treatment/exam table is needed. It is available in your choice of 20 colors of upholstery with color coordinated laminate doors and drawers on an easy-clean laminate base.

FEATURES

- Plush vinyl upholstery over 2" firm padding
- Full 27" width for patient comfort
- Air-spring adjustable backrest
- Coated, two position adjustable stirrups
- Padded, pull-out leg rest on steel track
- Optional pull-out basin
- Easy care all laminate base
- Color coordinated doors and top available
- Chip resistant edges
- 1 large side storage compartment with shelf
- Storage compartment is accessible from both sides (no right or left handed tables)
- 2 Euro style drawers with all steel slides
- 400 lbs. load capacity under normal use
- Built with plywood and laminate panels that are California CARB I compliant





Family Practice Exam Table

	Length	Height	Width
8870	58"/72" extended	31"	27"

- Full cabinet model table side storage access from either side
- 2 deep storage drawers on end
- Two position length adjustable stirrups
- Padded, pull-out leg rest
- Top sections are easily replaceable in the field
- Easy-clean laminate base
- Adjustable feet for leveling



Pull-out pad



2 deep drawers



Paper Roll Storage



2-Position Stirrups



Airlift Backrest



Both Sides Access




099 Optional Pull-Out Pan

ORDERING INFORMATION FOR MODEL # 8870

1. Select base frame color (page 4)
2. Choose a door and drawer color from the Color Selector on page 4
3. Pick a vinyl upholstery color from the Color Selector on page 4
4. See page 5 for drawer/door pull selection. If no selection is chosen Silver wire pulls will be used.

NOTE: If only base color is selected it will be used for the entire base.

TABLE OPTIONS



- *A 003 – Slide Out Shelf
Accessible from either side of table; 9" wide.
- B 004 – Drawer
Euro-style drawer with metal sides and laminate bottom and ends. Holds up to 75 lbs.
- *005 – Special Height
36" maximum, 21" minimum. When ordering, state desired height to top of pad on table.
- *006 – Special Length
78" maximum, 48" minimum length. When ordering state desired length (table top).
- *C 008 – Armboard
Adjustable in height and angle. Includes two mounting brackets.
- D 009 – Nosehole
Permits patients to lie face down without turning head. 9" long x 2" wide.
- *E 010 – Stirrups
Adjustable in height, angle, and length. Heels are coated; stirrup bars are chrome plated. Specify type of table when ordering.
- F 020 – Small Pillow
Dimensions: 14" x 12" x 3"
- 060 – Large Pillow
Dimensions: 14" x 21" x 3"
- G 030 – Paper Dispenser
Accommodates paper rolls from 18" to 21" wide.
- H 040 – Paper Cutter
Clear plastic with elastic straps and buckle
- I 077 – Safety Strap
2" x 92" Strap secures with quick fasten hook and loop fastener.
- J 088 – Flat Top Option
Flat table top option can not ship FedEx® Ground and MUST ship by common carrier with all other normal shipping charges
- *K 086 – Red Oak Option (Classic Series Only)
- *L 079 – Triple Bolted Legs Option (Classic Series Only)
- *M 018 – Removable Mat Top Option
2" thick removable mat is held in place with Velcro® strips (Only available with non adjustable backrest tops)

*option not available on Alpha Series tables.

Upholstery Colors	Base Laminate Colors
 3SB -Slate Blue	 1NT -* Natural
 3DT-Desert Tan	 1DC-* Dark Cherry
 3BK -Black	 1MP -Maple
 3PG -Purplegray	
 3VP -Viscaya Palm	
 3AS -Allspice	
 3BK -Black	
 3SP -Sapphire	
 3GP -Graphite	
 3AB -Aubergine	
 3AL -Alabaster	
 3NL -Neutral	

CLINTON VINYL UPHOLSTERY

Clinton uses only premium, woven knit-backed upholstery, that resists stretching and bunching associated with ordinary vinyl upholstery. Plus, Clinton upholstery is topcoated to resist stains, scuffs and abrasions. Standard flammability specifications meet CAL 117. For additional flammability requirements, call for quote.

ORDERING INFORMATION

- PICK A TOP COLOR**
Upholstery for table or laminate for cabinet
- CHOOSE A LAMINATE FRAME COLOR**
Legs, frame and table backs
- SELECT DOOR & DRAWER LAMINATE COLOR**
NOTE: If one color is selected, it will be used on the entire base
- SELECT DOOR & DRAWER PULLS**
NOTE: If no selection is noted on the order, Gray pulls will be used on solid colors, and Silver Wire pulls will be used on wood tones.



Color Coordinating Upholstery	Color Coordinating Laminates
 3GM -Gunmetal	 2SG -Slate Gray
 3DG -Dove Gray	 2GR -Gray
 3CM -Cream	 2FS -Fossil
 3CS -Clamshell	 2SD -Sand
 3WW -Wedgewood	 2AB -Aztec Blue
 3RB -Royal Blue	 2DB -Deep Blue
 3EM -Emerald	 2IG -Imperial Green
 3MB -Mulberry	 2WR -Wild Rose
 3BG -Burgundy	 2SGA -Sangria

Due to the limitations of the printing process, actual colors may vary. Actual upholstery and laminate samples are available.

Quick Ship Log—3 to 5 days
Quick ship items do have similar item quantity restrictions and also have some option limitations. Call for more information.

TREATMENT & BEDSIDE CABINETS

Available from 

FEATURES

Clinton Treatment Cabinets are designed with the needs of the professional in mind. Built tough with a wide range of sizes, models and a variety of storage configurations, Clinton treatment cabinets are the perfect complement to any office or clinic. Mix or match up to 12 colors of laminate, or coordinate with Clinton table upholstery.

ORDERING INFORMATION

When ordering cabinets you need to select top, front and side colors. If only one color is selected, that color will be used for the entire cabinet. One-piece ABS molded top color is always Gray.

- Easy-clean all laminate construction
- Tough chip-resistant edges
- Cabinets are assembled with steel fasteners, not glued wood dowels
- 3-way adjustable, 110° door hinges
- Self closing, Euro-style drawers with metal sides that hold up to 75 lbs.
- Recessed toe kick or dual wheel caster bases
- Adjustable drawer fronts
- Mix color choices for top, sides and fronts
- Select models are available with a 1-piece gray ABS molded top rimmed on all sides
- Built with laminated panels that are California CARB I compliant
- 200 lbs. load capacity

ABOUT PULL CHOICES

Clinton provides a choice of 2 contemporary styles of drawer and door pulls.

Durable plastic pulls in 4 complementary colors or a choice of brushed chrome or black wire style metal pull.

CONTEMPORARY AND WIRE PULLS



Gray 5G



Silver Wire 5C



Slate 5SG



Black Wire 5B



Sand 5S



Fossil 5F

If no pulls are specified, Grey pulls will be used on solid color fronts and Silver Wire pulls will be used on woodtones.

