

Datex-Ohmeda S5

Anesthesia Monitor

Equipment for the way *you* operate



 Certified Refurbished

All-inclusive monitor that is perfect for challenging anesthesia environments

With the Datex-Ohmeda S5 Anesthesia Monitor, you get a flexible range of monitoring parameters. The S5 also features module frame options/ displays and data transfer, and is easy to configure.

► A Quick Look

- User interface dedicated to anesthesia care, quick access to parameters and alarm management
- Supports up to eight waveforms
- Features user-configurable display layout and colors
- Up to 24-hour graphical and numerical trends
- Minitrends, split screen and four user configurable graphical trend pages
- Extended arrhythmia analysis and comprehensive range of advanced clinical parameters
- Quick user guidance via intuitive help function
- Six user-configurable modes with pre-configured settings

Overview of module frames and displays

- Available with either five or eight module slots; as many as four additional module slots can be added to the eight-slot configuration
- Insert/remove plug-in parameter modules with no interruption to other monitoring
- Removable ComBar and Remote Controller with Com-Wheel
- Various mounting elements for the system components