

DRE Anetic-P

Mobile Surgery Table (Powered)

Equipment for the way *you* operate

The complete multi-purpose mobile surgery table

Powered patient trolley/operating table is a total solution for hospitals and surgery centers.

The DRE Anetic-P is more than a surgery table. It is a complete surgery solution that can move patients from the pre-operative area to the operating room and then to the recovery area – with no need to transfer the patient from a stretcher to a surgery table and back. The Anetic-P does it all by incorporating the mobility of a stretcher into a multi-purpose surgery table.

Features

- » Four remote-activated powered functions
- » Lateral tilt
- » Excellent stability
- » Full C-arm coverage
- » Superior surgical access
- » Traversing patient platform
- » 550-lb. (250-kg) weight limit



Pictured with optional armboards



See the Anetic-P in action!

Visit dremed.com/anetic for a video that shows smooth positioning and essential accessories for the mobile surgery table!



Accessories

Surgical Headrests

Armboards

Leg and Foot Positioners

plus much more

Stretcher maneuvering is incorporated into the complete surgery system



Stretcher Mobility

Cotsides fitted for patient comfort and safety.

Retractable steering mechanism engaged for maximum maneuverability.



Full C-arm coverage

DRE Anetic-P

Mobile Surgery Table (Powered)

Equipment for the way *you* operate

Gynecology features

Patient platform traversed to foot end for access.

Head section removed for anesthesia access.

Leg section removed for access.

Fitted with optional powered leg supports.



± 12° Lateral Tilt

Shoulder Arthroscopy features

Patient platform traversed to foot end for access.

Optional arthroscopy backrest fitted.



Ophthalmology

Patient platform traversed to head end for access.

Ophthalmic head section fitted.



Removable head section



Push button rotating cotside



Push button release and self-locking leg section



Removable ultra lightweight leg section



K8 pressure care mattress



Technical Specifications

Height (min)	23.8" (605 mm)
Height (max)	39.5" (1,005 mm)
Trolley length	82.7" (2,100 mm)
Side Rail Width	25.8" (655 mm)
Cotside Width	30.3" (770 mm)
Brake Width (Brakes off)	33" (835 mm)
Mattress Length	80.3" (2,040 mm)
Mattress Width	23.6" (600 mm)
Mattress Depth	3" (75 mm)
Base Length	47" (1,190 mm)
Trendelenburg	20°
Reverse Trendelenburg	12°
Lateral Tilt	+12°/-12°
Backrest Articulation	0 - 65°
Head Section Articulation	+30°/-30°
Head End Traverse	9" (230 mm)
Foot End Traverse	9" (230 mm)
C-arm Access	100%
Max. Patient Weight Limits	550 lbs (250 kg)
Trolley Weight	352 lbs (160 kg)
Lightweight Leg Section Weight	13.2 lbs (6 kg)

Specifications subject to change without notice.

Trolley Features

- » Powered Functions
- » Lateral Tilt
- » Head End Traverse
- » Foot End Traverse
- » K8 Pressure Care Mattress
- » Removable Ultra-Lightweight Leg Section
- » Removable and Rotating Cotside
- » Push Button Release/Self Locking Controls
- » 5th Wheel Steering Mechanism
- » Advanced Four-Caster Braking System
- » E Size Oxygen Cylinder Trough (Side Loading)
- » Transfusion Pole with Mounting Position
- » 150 mm Castors

Optional Features

- » Spare Transfusion Pole
- » Sculpted Mattress for Head Support
- » Ophthalmic Head Support
- » Removable Shoulder Arthroscopy Backrest
- » Padded Cotside Covers - Gray
- » Padded Cotside Covers - Child Print
- » Articulating Removable Leg Section
- » Removable Leg Section with Full Length Side Rail
- » Press Mattress Complete
- » Accessory Trolley: Complete storage solution
- » Accessory Trolley Spare Trough
- » Spare Powered Trolley Hand Control
- » Spare Powered Trolley Charging Cable
- » CD Side Cylinder Support Bracket