

DRE Integra SP II (MRI-Compatible)

Anesthesia Machine

Equipment for the way *you* operate

DRE's MRI-compatible anesthesia machine is designed to meet all your anesthesia requirements

- Three gas machine (O₂ • N₂O • Air)
- Is a smart alternative to refurbished machines
- Advanced easy-to-use features for maximum patient safety
- Holds two vaporizers
- Large work area with a writing table and footrest

Order now

to receive a **two-year manufacturer's warranty** and a **two-year parts warranty**

► A Quick Look

- MRI-compatible
- Ventilator option (not included): Nuffield MRI-compatible anesthesia ventilator
- Selectatec®-compatible vaporizer mounting system
- Holds up to two vaporizers
- Specifically designed for low flow
- Left or right handed layout
- Contemporary ergonomic design
- Advanced easy-to-use features for maximum patient safety
- Mechanical anti-hypoxic device and air/N₂O interlock
- Only two preventive maintenances required each year

Physical Dimensions

- Size (h x w x d): 54" (137 cm) x 26" (66 cm) x 26" (66 cm)
- Weight (with all options): 200 lbs. (75 kg)

Power

- Power Requirements: 110/220 VAC
- Auxiliary Electric Outlets: (4) 13 amp

Display

- Top Shelf: 24" (61 cm) x 16" (40 cm)
- Work Surface: 24" (61 cm) x 14.5" (37 cm)
- Drawers: 6" (15 cm) x 15" (38 cm) x 17" (43 cm), maximum of three
- Writing Tablet: A4 size



Available with the Nuffield Anesthesia Ventilator

- MRI-compatible, pneumatically driven time-cycled ventilator
- Pre-set volume and flow rate for adult or pediatric patients
- Powered by medical air or oxygen at 340 kPa - 410 kPa

Ref. Code: U-8009

1800 Williamson Ct. • Louisville, Kentucky 40223 USA
800-477-2006 • 502-244-4444 • FAX: 502-244-0369

www.dremed.com

