

DRE Integra VS02

Portable Anesthesia Machine

Equipment for the way *you* operate.



the Compact and **Portable**

three-gas anesthesia system

The fully functional DRE Integra VS02 is compact enough for the smallest surgery or imaging suite. And at only 27 pounds*, moving the VS02 is a synch.



◀ Available with an optional king-bag

► A Quick Look

- Small and compact unit ideal for office based anesthesia
- Three gas machine (O₂, N₂O, and Air)
- Mechanical Anti-Hypoxic Device and Air/N₂O interlock
- Ventilator (options): DRE AV800, Ohmeda 7000 or 7800
- Selectatec®-compatible vaporizer mounting system
- Holds one vaporizer
- Specifically designed to accommodate low flow
- Auxiliary O₂ flowmeter port
- Recessed flush valve
- Optional carrying case
- Two-year manufacturer's warranty

On the move ►

- You can select one of two optional carrying cases:
1. With wheels
 2. Without wheels

Includes
a **2** year
manufacturer's
warranty

