

DRE LI-60

Tumescent Infusion Pump

Equipment for the way *you* operate

High-powered tumescent infusion pump designed for virtually effortless infiltration



The high-powered DRE LI-60 Tumescent Infusion Pump is designed for virtually effortless infiltration. It is capable of infusion flow rates of 0-850 ml/min and has an easy-load hub.

Features:

- Two infusion needles
- Handle with on/off control
- Infusion tubing (case of five)
- Pneumatic foot switch

Also available: 220 volt infusion pump

Attachments, Accessories and Supplies

Infusion Handle with On/Off Control

The new look of infiltration with on/off control. With an all stainless steel head control, the new design is for durability and a nice look. No tools are needed to disassemble this hand piece for cleaning.

Infusion Handle with On/Off Control Luer Lock Connect

The on/off valve workings have been improved by sealing systems that will better withstand the high autoclave temperatures. This improved handle allows for fingertip control during the infusion process.

Infusion Handle with Stop Cock Control

The level handle shows flow direction for easy visual checks, quarter turn shot-off. The valve is made of brass and chrome-plated finish. The stop cock control can be rotated to any direction by releasing the handle cap and tight the cap back after you choose what direction you want.

Luer Lock Infusion Needles

- For rapid infusion of tumescent fluid
- Standard diameters: 12, 14, 16 gauge
- Standard lengths: 15, 25, 30 cm

Garden Spray Infiltration Needle:

This cannula is used specifically for the infusion of tumescent solution. Configuration of the holes allow for even distribution of the fluid.

Super Spray Infiltration Needle:

The super spray is designed to maximize liposuction efficiency, with multiple round holes arranged in a helical pattern along the distal end of the cannula.

STD. Single Spike Infusion Tubing

- With silicone sleeve for infusion pumps
- Increased flow rate
- 188" length
- Available in a case of 10

Single Spike, Large Bore Tubing

- For manual and pressure cuff systems.
- Without silicone sleeve, 120" length.
- Available in a case of 10

Pressure Cuff for Infusion

- 1000 ml durable unit
- Easy to clean, cost effective and fast fluid level inspection
- 0-245 ml/min, constant pressure of 200 psi