

## Premium features and a sleek design meet affordability.

The **DRE Lucerne 360** provides intuitive positioning and includes a kidney bridge.

With the DRE Lucerne 360, you get premium features for as much as 30% less than the cost of competing surgery tables. Crafted for quality, the Lucerne 360 provides intuitive movement and stability, setting new standards in patient positioning. It features a C-Arm compatible table top with a 500-lb. (225-kg) lift capacity (also available: 1,000-lb. configuration).

The Lucerne 360 has many advanced positioning options, such as kidney bridge positioning, 360° horizontal rotation (180° each direction) and a horizontal slide motion that helps create a greater field of view for C-arm use.



**Rotates a total of 360°**  
for precise patient positioning



### Features

- » C-arm interface capabilities to meet imaging requirements
- » Remote control provides high operation safety and precise patient positioning with a safety lock
- » Silent electro-hydraulic movements
- » Mobile anti-static castors with a floor lock
- » Easy access to electro-hydraulic system for service
- » Kidney bridge and flex positioning
- » 500-lb. (225 kg) lift capacity – 1000-lb. (453-kg) option available

**Technical Specifications**

Bed size	2,050 mm x 500 mm 80.7" x 19.6"
Height	670 mm - 1,100 mm 26.3" - 43.3"
Reverse Trendelenburg	anterversion $\geq 30^\circ$ retroversion $\geq 30^\circ$
Title Left/Right	left oblique $\geq 30^\circ$ right oblique $\geq 30^\circ$
Chest Plate Adjustment	N/A
Back Plate Adjustment	+75° ~ -20°
Head Plate Adjustment	$\geq 40^\circ / \geq 90^\circ$
Leg Plate Adjustment	$\geq 90^\circ$
Leg Plate Angle	0 ~ 90°
Rotation	0 ~ 180° each direction
Longitudinal Shift	N/A
Wireless Remote Control and Effective Distance	N/A
Wired Remote Control	Standard
Four-level Position Reset Function	N/A
Emergency Power Supply	N/A

**Standard Accessories**

- » Anesthesia Frame (1)
- » Shoulder Supporter (2)
- » Arm Rest (2)
- » Leg-Crutches (2)
- » Waist Supporters (2)
- » Square Clamp II (2)
- » Square Clamp I (5)

**Optional Features**

- » Battery backup
- » Auto-reset / level function
- » Wireless remote control

