



## Cost-effective and manually-operated

Small details can bring big results, especially in today's surgical environment. DRE Marseille tables feature manually-operated functions for optimal performance and stainless steel columns that provide stability and a sleek look for any surgical suite.

This state of the art universal operating table combines the latest technology with the most affordable pricing in the market. The sturdy stainless steel construction and versatility allows the surgeon to perform many types of surgical procedures such as gynecological, obstetrics, ophthalmology, orthopedic, ENT, urology, abdominal and chest surgery.

The DRE Marseille is elevated by a hydraulic pump system operated with a foot pedal, making it extremely convenient and reliable.

All Trendelenburg functions and lateral movements are controlled by a mechanical operating system, which can be manually manipulated with handles conveniently located on the table's midsection.

The design and versatility of this universal operating table make a unique difference in the surgical suite.

## Key Features

- Manually operated functions
- Stainless steel construction
- Removable head and leg sections

## Specifications

- Tabletop length and width: 83" x 19" (210 x 48 cm)
- Elevation: 31 - 41" (80 - 104 cm)
- Trend / Reverse Trend:  $\pm 25^\circ$
- Lateral Tilt:  $\pm 18^\circ$
- Leg section adjustment:  $\pm 90^\circ$
- Back rest adjustment:  $\pm 75^\circ / -10^\circ$
- Head section adjustment:  $\pm 45^\circ / -50^\circ$
- Maximum kidney break position: 4" (10 cm)
- Maximum weight capacity: 381 lbs (173 kg)

## Standard Accessories

- Head rest: 1
- Anesthesia screen: 1
- Shoulder supports: 2
- Arm boards: 2
- Lateral supports: 2
- Stirrups: 2