

DRE Marseille

Manual Hydraulic Surgical Table

Equipment for the way *you* operate.



The manual, economical alternative to automatic tables

- Easy to adjust to Lateral, Trendelenburg and Reverse Trendelenburg positions
- Features a foot pedal for oil-hydraulic up and down positioning
- Dependable manual design makes the Marseille easy to maintain
- Standard accessories include a screen frame, head-rest, shoulder support, armrest, body support, knee crutches, and foot support
- Optional accessories include a restraint strap, hand grips, u-shaped headrest, raised armrest, stainless steel funnel and waste receptacle and corrugated tubing, and kidney elevator



The manual Marseille is easy to adjust to Trendelenburg, Reverse Trendelenburg and lateral positions

► Specifications

- Length of tabletop: 195 cm
- Width of tabletop: 55 cm
- Tabletop height (lowest/highest): 73 cm - 103 cm
- Lateral tilt (left/right): 25°/25°
- Trendelenburg/Reverse trendelenburg: 35°/35°
- Back section up/down tilt angle: 70°/30°
- Seat section up/down tilt angle: 70°/30°
- Leg section up/down/split tilt angle: 0°/90°/180°
- Head section up/down tilt angle: 25°/90°
- Weight: Approximately 170 kg

