

## DRE Millennium 3

Patient Stretcher

Equipment for the way *you* operate

The DRE Millennium 3  
is both **DURABLE**  
and **AFFORDABLE**



The DRE Millennium 3 is the durable stretcher alternative for budget-conscious facilities.

- Supports up to 550 pounds on a four-inch foam mattress
- Adjusts to Trendelenburg and Reverse-Trendelenburg positions
- Infusion pole sockets at each corner eliminate the need for additional equipment in small spaces and ensure a secure transport

### ► A Quick Look

- Twin pedal Trendelenburg
- Height adjustment clearly marked
- Pneumatic back adjustment
- Central locking casters
- Brake/steer pedals at each corner
- Fold down side rails
- Four infusion pole sockets
- Oxygen cylinder storage built-in
- 550-pound capacity
- 4" foam mattress
- Height adjustments from 23.5" to 35"
- Surface dimensions: 75" x 26" (l x w)
- Factory new with a 1-year warranty