

Strength and dexterity with a full-length radiolucent top.
Sleek, easy-to-clean stainless steel construction.

- » Dependable power with a 180° rotating top – provides reliable C-arm access
- » Completely radiolucent top, and superior imaging from an offset cylindrical column
- » Ample weight capacity – supports as much as 750 lbs.
- » Smooth steering castors provide enhanced mobility
- » Additional features and options: Battery back-up, brake system with auto-leveling, head and/or leg section (each removable), and optional traction bar



180° Rotating Top
for precise patient positioning



Specifications

Weight capacity

Lift	750 lbs.
Articulation	750 lbs.

Table top

Length	78"
Width	19.75"
Max. height	43"
Min. height	29"

Positioning

Lateral tilt	30°
Trendelenburg	30°
Back section	+90°, -40°
Head section	+30°, -90°
Leg section	+0°, -105°
Table rotation	180°
Self-leveling brakes	Powered



Easy to position by hand control or optional foot control