

Hill-Rom Advanta

For safety, quality and comfort.

Intelligent, well designed, clinically-effective and easy to use.

Refurbishing process includes

- Professionally refurbished medical beds to manufacturer's specifications
- JCAHO required documentation
- New paint & labels
- Fully function tested
- Reconditioned PC boards and side rail switches
- Reconditioned side rails
- Reconditioned head/foot boards
- Fresh hydraulic fluid
- Rebuilt casters with new components and wheels
- 1 year warranty (Domestic only)

Specifications

- Side communications (Nurse call, lighting, entertainment, volume and channel controls)
- Bed exit
- Zone air mattress
- Tente casters
- Scale
- Bumpers
- Automatic night light
- Instant CPR release
- Height in high position: 30"
- Height in low position: 17"
- Sleep surface width: 36"
- Sleep surface length: 82"
- Overall length (w/out bumpers): 36" x 91"

Options

- Service Contracts for preventative maintenance
- 220-volt transformers

