

MiniMax™ T3



NEVER COMPROMISE
START WITH YOUR VISION

Quartz-halogen light with matchless versatility that's compact and powerful



Remarkably illuminates
any headlight

MiniMax™ T3 is a powerful and compact 200-watt quartz-halogen light source that remarkably illuminates any headlight.

Offered with a three-port rotating turret – choose from ACMI, Wolf, Storz, Olympus and BFW – MiniMax™ T3 is outstanding for office use with numerous lighted instruments. Our benchmark manual intensity control quickly adjusts illumination and flips a filter in position for whiter light on demand. A cost-effective Lamp Life Saver Setting will double use time. Engineered with a thermostatically controlled fan, MiniMax™ T3 remains cool during and after use. An optional clamp-on tray for adaptability to your workstation or chair is available.



1-800-477-2006 • www.dremed.com



Our FOS 5T headlight system combines comfort and brilliant illumination with the High Lite Featherweight headlight and the compact, versatile and powerful MiniMax™ T3 light source.

Specifications

Product Codes	MID 2050 MiniMax T3
	FOS 5T MiniMax T3 Light Source With High Lite Headlight and 7' Cable
	MID 3210 Replacement Lamp
	MID 2010 Clamp-On Tray
	MID 2009 Mobile Stand
	MID 4005 Carrying Case for Light Source and Headlight
Port Options (Choose 3 When Ordering)	
	MID 2060 BFW
	MID 2061 ACMI
	MID 2062 Wolf
	MID 2063 Olympus
	MID 2064 Storz

Lamp	200-Watt Quartz-Halogen Color Temp. 3,600° Kelvin 4,000° Kelvin in White Filter Mode 40-Hour Lamp Life 80-Hour Lamp Life in Lamp Saver Setting
Dimensions	7.5" L x 5.5" W x 4.3" H (190 mm x 139.7 mm x 109.2 mm)
Weight	6.75 lbs. (3 kg.)
Power	100-240V - 50/60 Hz
Compliance	ETL Certified to UL 2601 CETL Certified to CAN/USA C22.2 No. 601.1 CE Certified

