



SPECIFICATIONS:

- weight: approximately 400 lbs
- max. dimensions: 40 x 68 x 25 in (l x h x w)
- top shelf: 32 x 15 in (l x w)
- table top: 271 in sq

Drager Narkomed 2C Anesthesia Machine

“The reassurance of a Drager without the expense of new...”

DRE offers the NAD Narkomed 2C refurbished to each customer’s specifications. This model holds 3 vaporizers for multiple agents and has oxygen, volume & pressure monitoring on a removable flat plate display. The manufacturer offers support of this model. Drawer pack configurations may vary from picture.

- O2 and N2O standard
- air optional
- ORC (Hypoxic Guard) standard
- AVE ventilator with ascending bellows standard
- holds 3 vaps
- compatible with 19.1 or 19.3 vaps, optional
- oxygen, volume & pressure monitoring standard
- removable flat plate display standard
- CO2 absorber standard
- E-Cylinder and central gas connections standard
- drawer pack configurations may vary from pic

