



Penlon Sevoflurane VAPORIZER

“...reflects the unique design and construction of the vaporizer, which gives greatly enhanced responsiveness to changes in gas flow during clinical procedures.”

The Sigma Delta is designed as a light weight service free* vaporizer for continuous flow anesthesia. It delivers accurate concentrations under varying conditions of flow rate and temperature, particularly at low flows.

Shown above as the Sevoflurane version. Halothane, Enflurane, Isoflurane compatible versions are also available.

This Penlon Vaporizer combines the latest patented technology and research to offer one of the highest quality options available.

| Vaporizer Type | Halothane 4% 5% 8% Filter Type | Enflurane 5% 7% Filter Type | Isoflurane 5% Filter Type | Sevoflurane 8% Filter Type |
|---------------------------------------|--------------------------------------|-----------------------------------|---------------------------------|----------------------------------|
| Selectatec® Compatible with Interlock | ✓ | ✓ | ✓ | ✓ |
| Drager Plug-In® Interlock Compatible | ✓ | ✓ | ✓ | ✓ |
| Cagemount Taper 23 mm | ✓ | ✓ | ✓ | ✓ |
| North American Drager | ✓ | ✓ | ✓ | ✓ |
| Selectatec® Compatible | ✓ | ✓ | ✓ | ✓ |

| | | | |
|------------------|--------------------------------------|-----------------------|---------------------------------------|
| | | | |
| Cagemount® Taper | Drager Plug-In® Interlock Compatible | North American Drager | Selectatec® Compatible with Interlock |

Technical Specification

Weight: Approximate weight 5 kg
Capacity: Volume at MAX mark: 250 mL, nominal Volume at MIN mark: 35 +/- 10 mL
Note: After draining, approximately 60 +/- 10 mL of liquid is retained by the wick
Flow Range: Operating flow range: 0.2 to 15 liters/min
Temperature Range: Operating temperature range: 15 to 39°C (59 to 95°F)
Cagemount Dimensions (mm): 133 wide x 158 deep x 219 high
Selectatec Compatible with Interlock Dimensions (mm): 120 wide x 150 deep x 242 high
Drager Plug-In Compatible Dimensions (mm): 100 wide x 158 deep x 242 high

Typical Performance Graph

PRODUCT OPTIONS

Agents: Halothane, Enflurane, Isoflurane, Sevoflurane

Mounting: Selectatec®, Drager® Plug-In, Cagemount

Filler: Keyed Filler, Quik Fil®, or Pour fill