



Welch Allyn Vital Signs 300 Patient Monitor

“The Welch Allyn Vital Signs Monitor 300 Series is an affordable vital signs monitor that’s both flexible and easy to use.”

Like the original Welch Allyn Vital Signs Monitor, the 300 Series helps enhance your productivity by providing automated vital signs before, during and after procedures. The simple, push-button menu operation and bright display allows you to quickly learn and easily obtain the important information you need, so you spend more time with your patients and less time with your equipment.

PRODUCT FEATURES

The Welch Allyn Vital Signs Monitor 300 Series begins as a non-invasive blood pressure (NIBP) monitor. Even in its most basic configuration, it delivers the features you need.

Multiple Patient Modes – Adult, Pediatric and Neonatal operating modes

Automatic Blood Pressure Mode – choose from 13 measurement intervals

Programmable Alarms – for high/low systolic, diastolic and pulse rate

LCD Display – Clock window and menu navigation

LED Display – Bright, clear viewing of vital sign readings

Built-In Memory – save up to 99 sets of measurement data

PRODUCT BENEFITS

- Separate patient modes with individual settings for adult, pediatric and neonatal support
- Fast, accurate temperature readings with industry leading Welch Allyn SureTemp® Plus thermometry
- Simple, straightforward operation and intuitive menu navigation with our innovative monitor design
- The Welch Allyn Vital Signs Monitor has a very short learning curve. Simple push-button programming and a large, digital display make learning and using the monitor easy and hassle-free.