

Amsco 3080

Surgical Table

The durable, flexible choice for your operating room

The Amsco™ 3080 features a radiolucent top and can be adjusted to Trendelenburg, back, side, foot and flex positions. It comes standard with a remote control for easy patient positioning.

A QUICK LOOK

- › Surgical table with radiolucent top
- › One of the most durable, flexible and reliable surgical tables available
- › Brings more efficiency to the operating room
- › Fluoroscopic capabilities
- › Simple patient positioning
- › Return to level with a single command

FEATURES

- › Remote control
- › Positions include: Trendelenburg, back (raise or lower), side (lateral tilt), foot (raise or lower), flex ("v" or inverted "v")
- › C-ARM Compatible
- › Return to level and special positioning models available
- › Up to 500 lb capacity
- › Refurbished with warranty



 Certified Refurbished

SPECIFICATIONS



Load Limit:
Up to 500 lbs



Maximum Height:
44"

Minimum Height:
27"



Trendelenburg:
25°

Reverse Trendelenburg:
25°

Lateral tilt:
18°

742-2018-06-26