

The Anyview 500R C-Arm offers vivid Images, innovative technology and unprecedented value. With intelligent low dose reduction, automatic brightness control and advanced pulse technology the new Anyview 500R is a cost-effective alternative to popular C-ARMS. Save money without sacrificing technology when you choose the Anyview 500R.



Anyview 500R Mobile C-Arm


Excellent radiographic imaging
at an economical price

The Anyview 500R Surgical Mobile C-Arm combines innovation and ease of use, making it the ideal option for a variety of applications.

- ✓ 9 or 12-inch High Resolution Image Intensifier
- ✓ Reduce radiation doses with pulse frequencies of 8, 15 or 30.
- ✓ Produces up to 30 frames per second in any imaging mode
- ✓ Kv and mA settings automatically optimize to the appropriate setting for each body part.
- ✓ Dual control panels allow clinicians to operate the C-arm from both sides of the unit.
- ✓ Intelligent Dose Reduction automatically reduces dosage without sacrificing image quality
- ✓ Tethered footswitch allows hands free imaging and image capture
- ✓ Dap software easily stores patient dose information and images
- ✓ DICOM 3.0 and PACS capable
- ✓ Stored Images on hard drive can be emailed, transferred to USB/DVD or local printer/network

C-Arm Comparison

Why Spend More?



	Anyview 500R	OEC 9800	OEC 9900 Elite	Philips BV Pulsera
X-ray tube anode	Rotating	Rotating	Rotating	Rotating
Radiographic mode	✓	✓	✓	✓
Fluoroscopic mode	✓	✓	✓	✓
17 kW @ 100 kVp	✓	✓		
Pulsed fluoroscopy	✓	✓	✓	✓
Pulses per second	8, 15, 30	1, 2, 4, 8	1, 2, 4, 8, 15	30
ABS control	✓	✓	✓	BodySmart
Snapshot mode	✓	✓	✓	Standard
Cine	30 fps	30 fps	30 fps	Standard
Reverse position	✓	✓	✓	✓
Warranty	2 year	1 year	1 year	1 year
C-Arm Measurements				
Free space	31"	31"	31"	30"
Depth	27"	26"	26"	24"
Orbital rotation	135°	148°	115°	135°
Horizontal travel	8"	8"	8"	8"
Vertical travel	19.7"	18"	18"	19.3"
Panning motion	12.5°	10°	10°	10°
Pivot rotation	360°	360°	360°	360°
Dimensions	70" x 32" x 71"	69" x 33" x 79"	69" x 33" x 79"	76" x 32" x 70"
Weight	552 lbs	720 lbs	606 lbs	683 lbs