

DRE ASG-200

Electrosurgical Unit

Pay for only the power you need

The DRE ASG-200 is a multipurpose electrosurgical generator for use in the modern operating room and surgi-center. It features both monopolar and bipolar functions to satisfy all of your surgical demands with safety, flexibility, reliability, and convenience.

FEATURES

- › 200 watts maximum power output
- › Monopolar and bipolar capable
- › 10 blend settings
- › Coagulation and fulguration modes
- › Pinpoint and spray coagulation
- › Split or solid return electrodes
- › Delivers consistent, repeatable power into varying load impedances
- › Defibrillator-proof type CF equipment
- › Footswitch indication on face of unit
- › Dual voltage 500/60 Hz
- › Does not require calibration
- › Two-year manufacturer's warranty



SPECIFICATIONS



Weight:
< 19 lbs (< 8.8 kg)



Height:
6.0 inches (15.3 cm)

Width:
12.3 inches (31.1 cm)

Depth:
16.3 inches (41.3 cm)



Cut:
200 W @ 300 Ω

Blend:
200 W @ 300 Ω

Pinpoint:
120 W @ 500 Ω

Spray:
80 W @ 500 Ω

Bipolar:
80 W @ 150 Ω

Duty Cycle:
10 S / 30 S

Output Frequency:
490 kHz ± 5

Input Current:
4.5 A ~

Line Voltage:
100 - 240 V ~

Line Frequency:
50-60 Hz

137-2018-06-27