

# DRE Avanti NXT3

## Infusion Pump

### Provides comprehensive infusion care for any medical facility

The Avanti NXT3 has an easy user interface allowing you to set all functions with just six buttons. It includes an extensive drug library as well as automatic and manual bolus.

#### FEATURES

- Avanti NXT3 Infusion Pump features varying modes of operation:
  - Vol/rate mode
  - Time mode
  - Body weight mode
  - Optional: Drop rate mode
- Onboard alarms indicate: infusion completion, low battery, air bubbles, downstream occlusion, no AC power, and more.
- Unique design allows for a more streamlined appearance and workflow.
- Pump design is compatible with straight line tubing sets from any manufacturer.
- Anti-bolus technology for increased safety after sudden occlusion release.
- Save and calibrate as many as three separate line sets.



#### SPECIFICATIONS



**Weight:**  
1.6 kg (3.5 lbs)



**Height:**  
193 mm (7.6 in)

**Width:**  
146 mm (5.74 in)

**Depth:**  
95 mm (3.74 in)

**Monochrome LCD Display:**  
3 in, 240 x 128 pixels



**AC Power Supply:**  
100-240 V; 50/60 Hz

**External DC Power Supply:**  
DC 10V-16V



**Mechanism:**  
Linear Peristaltic Movement

**Preset volume:**  
1-9999 ml

**Accumulated Infusion Volume:**  
1-9999 ml

**Increment:**  
0.1 ml/h

**Accuracy:**  
± 5%

2079-2018-10-25