

DRE Delphi CF

C-Arm Table

A 5-movement C-arm table with floating top

The DRE Delphi CF is an affordable C-arm tables with high-end, easy-to-use features.

FEATURES

- › Easy-to-use, five-movement C-arm table with high-speed actuators.
- › Rapid positioning via five table motions: height, lateral tilt, Trendelenburg, longitudinal travel, lateral travel.
- › Carbon fiber table allows clinicians to get high-quality, clear images with low X-ray doses.
- › Micromotors ensure quick, smooth, and balanced table movements.
- › Comes standard with a variety of accessories.
- › Maximum patient weight of 440 pounds. A 600lb capacity version is available. Contact your DRE representative for details.
- › Comes standard with onboard control panel and handheld remote control.
- › High-quality casters and braking system allow clinicians to easily move the table.
- › Compatible with a variety of popular C-arm systems from industry-leading brands.

STANDARD ACCESSORIES

- › Mattress - 1 piece
- › Onboard Controller - 1 set
- › Remote Control - 1 set
- › Lateral support - 1 pair
- › IV Pole
- › Fixed belt - 1 pair
- › Arm boards - 1 pair
- › Anesthesia screen
- › Removable rail - 2 pieces



1515-2018-08-22

SPECIFICATIONS



Maximum Patient Weight:
440 lb (200 kg)
Optional: 600 lb (272 kg)



Tabletop:
96.75"

Width of Tabletop:
20.9"

Imaging Area:
72.5" (head section)
32.75" (foot section)

Trendelenburg:
0-27°

Reverse Trendelenburg:
0-27°

Lateral Travel:
8"

Lateral Roll:
0-27°

Longitudinal Travel:
31.5"

Vertical Travel:
16.75"



Power:
120 V / 50 Hz / 350 VA
(available in 220 V)

OPTIONAL ACCESSORIES

- › Lead drape