

DRE HydroTherm SG

Autoclave

Operates easily, safely, and reliably

The DRE HydroTherm SG is a state of the art autoclave, designed for clinical and laboratory use. This autoclave includes an easy to use interface, intelligent control systems that prevent errors, and five pre-set options for convenience.

FEATURES

- › Five pre-set options for convenience:
 - › 250°F/solid (250°F/unwrapped)
 - › 250°F/wrapped (250°F/wrapped)
 - › 273°F/solid (273°F/unwrapped)
 - › 273°F/wrapped (273°F/wrapped)
 - › Cotton
- › Screen display of parameters and conditions during processing
- › Adjustable sterilization times
- › Equipped with a water quality check system to prevent dirty water being used
- › Waste water tank indicator light detects when water levels are too high/too low
- › Successive sterilization process option keeps the autoclave warm until it is opened, allowing for reduced warming time between sterilizations
- › If an error occurs, the machine will stop running, pressure will be released, and an alarm will go off to prevent damage
- › Key for vacuum test program



SPECIFICATIONS



Height:
16.5"

Width:
19.75"

Length:
25.5"

Chamber Dimensions:
9.75" x 17"



Frequency:
60Hz

Voltage:
one-phase
110V±11V



Volume:
Holds 22 liters

Components

- › 3 shelf rack
- › Equipment tray
- › Tray shelf
- › Tray holder
- › Drain hoses

2072-2018-10-09