

# DRE Meridian

## Power Exam Table

### Designed for patient comfort and convenience

The DRE Meridian Exam Table is an economical option for facilities seeking a simply designed, convenient exam table. The table comes with a pull-out foot rest and extendable stirrups, making it ideal for a variety of office-based medical exams. With electronic height adjustment, adjustable foot stirrups, and convenient paper roll holder, the Meridian Exam Table can be easily prepared between patients in any busy practice.

#### FEATURES

- › Ergonomically designed seat section and back section are designed for patient comfort.
- › Seat section features pelvic tilt for enhanced positioning capabilities.
- › Convenient and flexible operation with quiet motor.
- › Easy-to-clean, smooth surface.
- › Foot stirrups angle can be positioned for each patient: 15°, 30°, 45°
- › Extendable, padded foot rest.
- › Use convenient handle to manually adjust back and seat rest.
- › Electronic height adjustment via hand control.
- › Stainless steel basin slides out from the base of the table.
- › Convenient paper roll holder located under the back section.
- › Height adjustment range of 16 - 35.5 inches.



1562-2018-04-27

#### SPECIFICATIONS



**Maximum load:**  
550 lb



**Dimensions:**  
55 in ± 0.8 in x 28 in ± 0.8 in

**Minimum Height  
(without cushion):**  
16 in

**Maximum Height  
(without cushion):**  
35.5 in

**Cushion height:**  
3.5 in

**Slide length of foot rest:**  
11 in ± 0.8 in

**Slide length of pull-out basin:**  
17 in ± 0.8 in

**Table height adjustment:**  
19 in ± 0.4 in



**Back section angle:**  
0° ~ 70° ± 5°

**Seat section angle:**  
0° ~ 7° ± 1°

**Foot rest angle:**  
45° ± 5°