

DRE Millennium 7

Hospital Stretcher

Efficient patient stretcher with convenient, fifth-wheel steering

The DRE Millennium 7 is an industry-standard hydraulic patient stretcher designed for patient comfort and caregiver convenience. It's advanced steering system and ergonomically positioned pulling handles make it easy to maneuver. The durable build of the Millennium 7 ensures a safe and efficient stretcher that will meet the needs of any healthcare facility.

FEATURES

- › Integrative ABS plastic base
- › Convenient storage for patient belongings
- › Foot pedal controls on both sides for the operator
- › Dual pneumatic pistons aid in back section adjustment
- › Radio translucent backrest with X-ray cassette holder clip
- › Easy locking controls ensure patient and operator safety
- › Adjustable four-hook IV pole on foot and head ends of stretcher
- › Oxygen cylinder holder
- › Side rails fold down for easy patient access
- › Fifth wheel for easy maneuvering
- › Bumper casters on each corner to avoid collision damage
- › Large central locking casters on foot and head ends
- › Adjustable knee gatch for increased patient comfort



1961-2018-05-24

SPECIFICATIONS



Weight Capacity:
350 lbs (160 kg)



Total Length:
87"

Total Width:
30.5"

Platform Surface:
75.5" (L) x 24.5" (W)

**Height Adjustment
Range:**
22" - 35"



**Trendelenburg/
Reverse Trendelenburg:**
0 - 10°

**Back Section
Adjustments:**
0 - 80°

**Knee Gatch
Adjustments:**
0 - 35°