

DRE SL

Patient Exam Table

Designed for comfort, economy and overall functionality

This versatile table provides a simple solution for family practice, doctor's offices, clinics or any facility seeking a fully-functional, cost-effective exam table. The DRE SL Patient Exam Table features adjustable stirrups, plush vinyl upholstery. With chip-resistant edges and easily replaceable top padding, this table is built to last.

Features

- ✓ Plush vinyl upholstery over 2" firm padding
- ✓ Full 28.5-inch width for patient comfort
- ✓ Air-spring adjustable backrest
- ✓ Coated, two position adjustable stirrups
- ✓ Pull-out, laminate leg rest on steel track
- ✓ Compact design for economical shipping
- ✓ 400-pound weight capacity under normal use
- ✓ Easy care all laminate base
- ✓ Chip resistant edges
- ✓ Built with plywood and laminate panels that are California CARB I compliant
- ✓ Top sections are easily replaceable in the field



466-01-2017-02

SPECIFICATIONS



Length:
58in / 72in extended

Width: 28.5in

Height: 31in



Features a streamlined, straightforward design and comes in a compact package, allowing you to save on shipping costs

