

DRE Waveline ECO

Patient Monitor

User-friendly, high-resolution monitor offers superior performance

The cutting-edge innovation of the DRE Waveline ECO measures up to six parameters on a high-tech, low-power screen. With six easy-to-use interfaces to choose from, this advanced patient monitor is perfect for any clinical setting.

The DRE Waveline ECO is the top choice for health care professionals who demand high performance and precision, and is available through DRE Medical.

FEATURES

- › Designed for use with adults, children, and neonates.
- › Optional parameters: EtCO₂, Nellcor SPO₂, Suntech NIBP/IBP.
- › Versatile mounting options include: bedside hook, wall-mounted hanging bracket, mobile monitor stand.
- › Large storage function for NIBP measurement, alarm event, trend diagram/table review, etc.
- › Large font display to offer clear view of all vital signs.
- › Convenient touchscreen lock-out function.
- › TFT color touchscreen provides a brighter screen with better contrast and better visual effects, while using less power compared to LED screens.
- › Enhanced storage compared to other leading monitor brands:
 - › 2000 NIBP measurements
 - › 200 alarm events
 - › 48 hours waveform data
 - › 1000 hours trend diagram review
 - › 1000 hours trend table review
 - › 120 minutes dynamic trend

SPECIFICATIONS



Display:
12.1"

Resolution:
800 x 600



Battery:
3.5 hour
battery duration

Interface Options

- › All leads display interface
- › Trend coexistence interface
- › Large font Interface
- › NIBP review interface
- › Respiratory oxygenation map interface



1958-2018-06-28