



640

MIDMARK PEDIATRIC EXAMINATION TABLE





The Midmark 640 Pediatric Examination Table offers all the features and benefits you need to give your patients the best, most tender care possible.

With an upholstered end panel and beveled sides, infants can't roll off the Midmark 640. And, because the sides are flat with no sharp edges, you have peace of mind knowing your patients are safe from harm.

Bending to examine 25-30 children every day can put strain and stress on your back. The Midmark 640 eliminates that strain and stress with its exam surface height of 36". You get great access to the patient without the need to bend over and do your exam.

When you're ready to examine a patient, the Midmark 640 provides everything you need right at your fingertips.

There's no need to walk to another area for measuring or weighing, or to retrieve supplies. And there's no need for a separate scale that can take up valuable counter space.

Measure babies easily and accurately with the 640's infantometer. It pulls out easily and articulates out of the way when you are finished measuring.

For all your storage needs, the Midmark 640 provides one oversized drawer, as well as a bulk storage area that's perfect for keeping diapers and equipment out of the way. A hidden paper roll holder, with a paper tear strap at one end, is yet another convenience. For security, four leveling screws on the 640's base make it rock solid...and more importantly keep children from being able to crawl beneath the table.

Specifications

Working height:
36" (91.4 cm)
Table top length:
42" (106.7 cm)
Table top width:
20.75" (52.7 cm)

Paper roll holder

Upper portion of table accommodates one 14" x 3" (35.6 x 7.6 cm) diameter

Drawers

One large drawer:
33" x 15" x 4.125"
(83.8 x 38.1 x 10.5 cm)

Drawer is constructed of high impact polystyrene with rounded corners for easy cleaning and is supported by side-mounted ball bearing steel slides.

Bulk storage:
33.5" x 17" x 18.38"
(85 x 43.2 x 46.7 cm)

Options


Door/drawer locks



Built-in Scale

The built-in scale is easy to use. Simply turn it on, wait a second for it to calibrate, and it's ready to weigh. What could be easier? The unit comes with an AC adapter, or, if you prefer, it can be battery-operated. The unit even turns itself off, automatically.


Midmark Premium Solid Upholstery Colors



Midmark Designer Pattern Upholstery Colors



Ritter Upholstery Colors



Midmark is an ISO 13485 and ISO 9001 Certified Company. For more information or a demonstration, contact your Midmark dealer or Call: 1-800-MIDMARK Fax: 1-800-365-8631 Outside the U.S.A. Call: 1-937-526-3662 Fax: 1-937-526-8392 or visit our website at midmark.com

© 2008 Midmark Corporation
60 Vista Drive P.O. Box 286 Versailles, Ohio 45380-0286
Products subject to improvement changes without notice
Litho in U.S.A. 007-0674-00 Rev. A1 (1/08)



Because we care.