

Penlon Sigma Delta

Anesthetic Vaporizer

Offers internationally recognized performance and value

Penlon Sigma Delta vaporizers are available in a variety of agent and mounting options from DRE Medical. Their signature, long-lasting design offers clinicians peace of mind at an economical price point.

Features

- ✓ Anesthetic vaporizer offers excellent performance, especially with low flows.
- ✓ Internationally recognized, award-winning design.
- ✓ Multiple agent filler block options: Sevoflurane, Isoflurane, Halothane and Enflurane.
- ✓ Four back-bar mounting systems: Selectatec®, Dräger Plug-In®, North American Dräger, and Cagemount.
- ✓ Low body weight for easy use.
- ✓ Service-free vaporizer designed for years of regular use with without preventive maintenance services
- ✓ Key fill or pour fill models available — ask your DRE representative for details!



1065-01-2017-04

SPECIFICATIONS



Weight:
4.8kg approx



Cagemount Dimensions:
133 W x 158 D x 242 H (mm)

Selectatec® Compatible with Interlock Dimensions:
120 W x 190 D x 242 H (mm)

Dräger Plug-In® Compatible Dimensions:
100 W x 190 D x 242 H (mm)



Capacity Volume:
at MAX mark
250 mL nominal

Volume at MIN mark:
35 ±10 mL – Note after draining, approximately 60 ±10 mL of liquid is retained by the wick

Flow Range: Operating flow range 0.2 to 15 litres/min

Temperature Range:
Operating temperature range 15° to 35°C (58° to 95°F)

**Presumes ten year product life requiring no preventive maintenance service. It is recommended that a service is carried out at ten years. Halothane Vaporizers require a five year service.*

Selectatec® is a Datex Ohmeda Trademark, Plug in® is a Dräger Trademark