

SALTER LABS®



Salter-Style® Adult

Cannula

The Salter 1600Q (for "Quiet")

For oxygen flows up to 8 Lpm.

The 1600Q's larger reservoir virtually eliminates noise, creating a quieter, more soothing atmosphere.



Features

- Incorporates a larger reservoir
- Curved, tapered nasal prongs
- Three-channel, crush-resistant safety tubing
- Same unique anatomically-correct design as the Salter 1600

Benefits

- Virtually eliminates noise
- Better anatomical fit for long term use
- Uninterrupted flow of oxygen
- A more comfortable fit


Try the Salter 1600Q and HEAR the difference!



High flow oxygen patients now can enjoy quiet as well as smooth delivery

The Salter 1600Q Cannula

Clear, lightweight cannula with unique one piece anatomical design, eliminates the need for irritation-producing flap, result: maximum patient comfort. Over-the-ear style ensures secure positioning, allowing complete freedom of movement. Available with three-channel safety tubing or smooth bore tubing.



Ordering Information

Adult Nasal Cannula*

	Units (per case)	Safety Channel	Smooth Bore
Salter-Style® Cannula			
Clear cannula with extra-large facepiece for high flow situations, with 7' supply tube	50	1600Q	1650Q
Half case packaging of 1600Q cannula	25	1600Q HC	
Clear cannula with extra-large facepiece without supply tube	50	1606Q	



100 W. Sycamore Road - Arvin, CA 93203
 Phone: 661-854-3166
 1-800-421-0024 (CA only)
 1-800-235-4203 (U.S. & Canada)
 Fax: 661-854-3850
 Toll Free Fax: 1-800-628-4690 (U.S. & Canada)
 International Fax: 661-854-3850
 E-mail: salterlabs@us.salterlabs.com
 www.salterlabs.com



0482