



“Art is the expression of the profoundest thoughts in the simplest way”

---- EINSTEIN





5000 Hydraulic Patient Chair
(Shown with plush upholstery including trim)

- | | |
|---|------------|
| 5000 Hydraulic with Plush Upholstery and Trim | 5000-001T |
| 5000 Hydraulic with Sling Plush Upholstery and Trim | 5000-001S |
| 5000 Hydraulic with Sling Asepsis Upholstery | 5000-001SA |
| 5000 Hydraulic with Asepsis Upholstery | 5000-001 |
- (Plush upholstery available in Ultraleather only)*



Premier Orthodontic Patient Chair
(Shown with D5000 Ultra Doctor Stool and plush upholstery in Ultraleather)

- Premier Ortho Chair with Plush Upholstery and Trim
 - Premier Ortho Chair with Sling Plush Upholstery and Trim
 - Premier Ortho Chair with Sling Asepsis Upholstery
 - Premier Ortho Chair with Asepsis Upholstery
- (Plush upholstery available in Ultraleather only)*



2001 Electro-mechanical Patient Chair
(Shown with plush upholstery including trim)

- | | |
|--|------------|
| 2001 Electro-mechanical with Plush Upholstery and Trim | 2000-401T |
| 2001 Electro-mechanical with Sling Plush Upholstery and Trim | 2000-401S |
| 2001 Electro-mechanical with Sling Asepsis Upholstery | 2000-401SA |
| 2001 Electro-mechanical with Asepsis Upholstery | 2000-401 |
- (Plush upholstery available in Ultraleather only)*



2001E Economy Patient Chair
(Shown with optional Interceptor and Asepsis Upholstery)

2001E Economy Electro-mechanical with Asepsis Upholstery



cluding trim)

- Trim 2000-060T
- Upholstery and Trim 2000-060S
- Upholstery 2000-060SA
- 2000-060
- (Ultraleather only)




Patient Chair
(International Package)
(Upholstery)

Upholstery 2000-017



OS III Oral Surgery Chair
(Shown with optional Eight-way I. V. Arm)

OS III Oral Surgery Chair with Standard Upholstery 2000-080-2



OS VIII Oral Surgery Table
(shown with optional Chinner Headrest)

OS VIII Oral Surgery Table with Standard Upholstery 2000-088

ECONOMY WITHOUT COMPROMISE Craftsmanship of the highest order is a time honored tradition at Westar. The latest technology is combined with the skill of craftsmen to create hallmark products which are responsive, and assured to bring you affordable luxury. This makes Economy Without Compromise the tradition of choice.

5000 Hydraulic and 2001 Standard Features:

Wheelchair accessible double articulating headrest
 7 way foot control with two programmable preset positions and automatic return (exit)
 Metal construction with HiGlos powdercoat inside and out
 Quick Change Upholstery
 Synchronized seat and lumbar with back recline
 Thermal overload protection
 Low voltage
 Left/Right patient access with arms that swing up
 Locking 60 degree swivel
 Three year warranty
 Lifetime bushings throughout all pivot points

Additional features included with 5000 Hydraulic Patient Chair:

Hydraulic cantilever system
 Safety shut-off plate

Additional features included with 2001 Patient Chair:

Electro-mechanical cantilever system

Premier Ortho Patient Chair Standard Features:

New style Non-articulating headrest
 Base mounted Foot switch with manual operate and automatic return (exit)
 Pedestal base with 5" manual height adjustment
 Metal construction with HiGlos powdercoat inside and out
 Quick Change Upholstery
 Synchronized seat and lumbar with back recline
 Thermal overload protection
 Low voltage
 Three year warranty
 Lifetime bushings throughout all pivot points
 No arms

Optional features:

Fixed arms	2000-243
Movable Arms (<i>arms swing up</i>)	2000-244
Upholstered covered arm pads	7000-096-1

Optional Features for all chairs:

Ponytail headrest	7000-609
Non-articulating headrest (<i>Can be purchased in lieu of</i>)	7000-028
Double articulating headrest	7000-021
Touch-pad chair control with 5 programmable presets (<i>Upgrade for 2001E Economy chair required</i>)	
Touchpad for chair control	
5000: Chair mount	2000-230
Delivery mount	2000-233
2001/OS III: Chair mount	2000-231
Delivery mount	2000-234
Chair Adapter for 2" pole	2000-226
Internal umbilical	4000-054-1
Upholstered arms for non-plush upholstery	7000-096-1
Ultraleather upgrade - chairs (<i>not plush styling</i>)	1005
Ultraleather upgrade - stools and small parts	1006
Florentine Burl Trim	2000-229
Extended Three Year Warranty	1008
Child Seat Wedge	2000-246
Pedo upholstery (<i>Asepsis style Narrow back 6" shorter seat</i>)	7005-716

2001 Economy Standard Features:

Wheelchair accessible double articulating headrest
 Electro-mechanical cantilever system
 5 way foot control with auto-return
 Metal construction with HiGlos powdercoat inside and out
 Quick Change Upholstery
 Synchronized seat and lumbar with back recline
 Thermal overload protection
 Low voltage
 Fixed arms
 Three year warranty
 Lifetime bushings throughout all pivot points

OS III Oral Surgery Chair Standard Features:

Left/Right upholstered I. V. arms with one interchangeable pillow
 Wheelchair accessible double articulating headrest
 Adjustable / Removable arms
 7 way foot control with two programmable preset positions and automatic return (exit)
 Metal construction with HiGlos powdercoat inside and out
 Quick Change Upholstery
 Synchronized seat and lumbar with back recline
 Thermal overload protection
 Low voltage
 Three year warranty
 Lifetime bushings throughout all pivot points
 Electro-mechanical cantilever system
 Floor to seat: 13" to 31"

OS VIII Oral Surgery Table Standard Features:

Left/Right upholstered I. V. arms with one interchangeable pillow
 Wheelchair accessible double articulating headrest
 Flat when fully reclined
 Toe articulation
 Locking 60 degree swivel
 Narrow back upholstery
 Adjustable / Removable arms
 7 way foot control with two programmable preset positions and automatic return (exit)
 Metal construction with HiGlos powdercoat inside and out
 Quick Change Upholstery
 Synchronized toe with back recline
 Thermal overload protection
 Low voltage
 Three year warranty
 Lifetime bushings throughout all pivot points
 Electro-mechanical cantilever system
 Floor to seat: 20" to 35"
 Floor to Headrest: 43" to 58"
 Fully reclined (headrest to end of toe): 60"

Optional Features OS III and OS VIII:

Chinner headrest	7000-608
Eight way articulating I. V. arm board	2000-220
3" Height adjustment spacer	2000-119
Leg/Body Strap set	2000-223
Mobile Base	7000-371

Made in the U. S. A. 05/09