




Certified Refurbished

Easy-to-use infant incubator with a simple, easy-to-read display



- Provides stable temperature
- Mattress tray tilt mechanism
- Variable height adjustment
- Baby-temperature and air-temperature controls

Specifications

Dimensions: 40.4 mm x 119.4 mm x 59.4 mm

Weight: 140.6 kg

Power Supply: 240 V 50 Hz ac

Access

- Hand ports: 5
- Doors: 1 (front)
- Tubing ports: 6
- Tray withdrawal: Front
- Tray tilting mechanism: Pneumatic $\pm 9^\circ$

Controller

- Air-temperature: 20° C to 39° C (0.1° C)
- Baby-temperature: 34° C to 38° C (0.1° C)

Alarms

- Mains power failure
- High air temperature
- Skin temperature