



Certified Refurbished



## Flexible positioning and superb optical performance

Amsco Quantum SQ-240 surgical lights combine the flexibility and positioning ease of small diameter lightheads with the superb optical performance of large diameter lightheads. They are high-powered, easy to maneuver, and designed for most surgical procedures.

Configurations are available with a self-contained wall-mount intensity controller and optional sterilizable handles.

### ► Specifications

- Intensity levels: Up to 12,000 footcandles or 129,170 lux
- Cool illumination: Optical coatings remove 98% of all unwanted heat
- Color correction: Up to 4200° Kelvin
- Color rendering index: 92
- Continuous rotation: 360° at all critical joints
- Lateral illumination: 90°
- Continuous light: If a lamp fails, system automatically changes to spare lamp
- Lamp life: ~1,000 hours typical; 2,000 on Setting 3
- Depth of field: 26"
- Cylinder of light: 18"
- Head size: 24"
- Peak intensity: 18"