

Avante Torino EXL

Surgery Table

The Avante Torino EXL automatic operating table is an excellent choice for use in a variety of operations. A micro touch remote control operates the quiet electric motors, flexible adjustments of the head, back, and seat plates, and the built-in kidney bridge. The Torino EXL features long horizontal sliding, C-arm and x-ray access.

FEATURES

- › C-arm and X-ray interface capabilities to meet imaging requirements
- › Micro touch remote control with low noise movements
- › Flexible adjustments on head plate, back plate, seat plate and built-in kidney bridge
- › Quiet and reliable electric motors
- › Stainless steel design for excellent stability



512-2019-10-02

SPECIFICATIONS



Length/Width:
2070mm/550mm;
81.5in/21.7in

Elevation (down/up):
700mm/1000mm;
27.6in/39.4in

Kidney Bridge Up:
110mm; 4.3in

Horizontal Sliding:
300mm; 11.8in



Trendelenburg: 25°

Lateral Tilt (left, right): 15°

Head Plate (up/down):
45°/90°

Leg Plate (up/down/outward):
15°/90°/90°

Back Plate (up/down):
75°/20°



Power:
220V, 50 Hz, 1Kw

ACCESSORIES

- › Anesthesia screen
- › Shoulder support, 1 pair
- › Body support, 1 pair
- › Arm support, 1 pair
- › Leg support, 1 pair
- › Foot plate, 1 pair