



Certified Refurbished

Designed to meet the fluid and drug delivery requirements of today's changing medical environment

The AS40A Infusion Pump provides accurate continuous and intermittent infusion via intravenous (IV), intra-arterial (IA), epidural and subcutaneous routes of administration. It accepts standard disposable syringes from 1 mL to 60 mL in size. It also has a numeric keypad that simplifies programming and makes the pump easier to use.

Safety and effectiveness are reinforced by pre-programmable bolus operation, titration of a dose without interruption of fluid flow, and easy-to-understand alarm and alert messages.

The AS40A can run on its internal rechargeable battery pack and can also be operated while attached to a battery charger. It is supplied with a pole and built-in IV pole loop, and can be used as a table-top unit.

► A Quick Look

- Continuous or intermittent infusion
- Standard disposable syringes from 1 mL to 60 mL
- Numeric keypad
- Pre-programmable bolus operation
- Alarm and alert messages
- Internal rechargeable battery pack
- Multiple configurations
- Built-in IV pole loop



Tailor the AS40's key features to your requirements.

Configurable options include:

- Drug library feature
- Syringe manufacturer
- Automatic syringe size recognition
- Selectable infusion modes
- Maximum infusion rates
- Occlusion pressure sensitivity
- Keypad auto lock.

Selected options can be adjusted to meet new or different requirements.