



An essential resource in outpatient surgical settings

Administers medications that require multiple dosage adjustments

The Infus O.R. syringe infusion device aids in administering many intravenous agents given during anesthetic procedures. It conveniently delivers narcotics, muscle relaxants and vasoactive drugs via SMART LABELS—a system that displays appropriate infusion parameters and rates for specific medications, helping decrease the likelihood of calculation errors in fast-paced surgical environments.

► A Quick Look

- Modify dose without multiple reprogramming steps
- Adjust settings easily at any time during infusion
- Change dial settings with a single hand for an immediate response
- Save calculation time
- Minimize risk of calculation error
- Change infusion parameters/rates with SMART LABELS
 - SMART LABELS are magnetically coded/drug specific
 - Select from the available SMART LABEL library
 - Ideal for administering drugs requiring multiple dosage adjustments