

Perfect for in-office, minor skin-related procedures

Save your practice valuable time and money by using the Bovie DERM 101 10-watt high frequency monopolar desiccator to perform in-office, minor skin-related procedures. The DERM 101 comes ready to use with a 3-Button Power Control Handpiece to switch between power levels in 1/10th watt increments during procedures.

Features:

- Ideal for the in-office removal of skin tags, premalignant lesions, benign lesions
- Provides up to 10 watts of power
- 3-Button Handpiece performs desiccation and fulguration in 1/10th watt power increments
- On-board rotary dial may also be used to adjust power settings
- Bright LED display for clear readings anywhere in the procedure suite
- Unit comes ready to use, packaged with a selection of sterile and non-sterile tips

Standard Configuration:

- DERM 101 Unit
- 3-Button Power Control & Autoclavable Handpiece
- Sharp Dermal Tips, Non-Sterile (5 ea)
- Sharp Dermal Tips, Sterile (2 ea)
- Blunt Dermal Tips, Non-Sterile (5 ea)
- Blunt Dermal Tips, Sterile (2 ea)
- Derm-Elite™ Premium Blunt Dermal Tip, Non-Sterile (2 ea)
- Derm-Elite™ Premium Blunt Dermal Tip, Sterile (2 ea)
- Disposable Handpiece Sheath, Non-Sterile (2 ea)
- Wall Mount Kit (1 ea)
- Medical Grade Power Cord
- User Manual (CD Format)

Technical Specifications:

- Overall Dimensions: 2.7" (69mm) H x 7.2" (184mm) W x 8.2" (208mm) L
- Weight: Less than 3 lbs | 1.4 kg
- Voltage: 100 - 240 VAC ± 10%
- Frequency: 50/60Hz
- Current: 320mA Max
- Power Consumption: 28 VA

