

BARIATRIC SERIES



Bariatric Full Cabinet Treatment Table

	Length	Height	Width
6160	78"	31"	33"

- Plywood top secured in place with tee nut fastening system
- Steel reinforced stress points
- Cabinet constructed with steel fasteners
- Furniture-grade laminated plywood cabinet frame
- 3 storage areas, each with an adjustable shelf
- 6 doors with concealed, self-closing Euro-style hinges
- 600 lbs. capacity under normal use



Bariatric Treatment Table

	Length	Height	Width
6190	78"	31"	30"

- Plywood top and solid hardwood legs
- 6-leg design with triple bolted corner legs
- Built-in middle leg levelers
- Steel reinforced frame
- Laminate full length shelf with H-braces
- Pillow included
- 600 lbs. capacity under normal use



Bariatric Treatment Table/Low Height

	Length	Height	Width
6196	78"	26"	30"

- Plywood top and solid hardwood legs
- 6-leg design with triple bolted corner legs
- Built-in middle leg levelers
- Steel reinforced frame
- Laminate full length shelf with H-braces
- Pillow included
- Lower height makes patient access easier
- 600 lbs. capacity under normal use




Bariatric Treatment Table

	Length	Height	Width
6198	72"	31"	30"

- Plywood top and welded steel frame
- 6-leg design with bolted-through assembly
- Leg levelers
- Laminate full length shelf with H-braces
- 600 lbs. capacity under normal use

TABLE OPTIONS



- *A 003 – Slide Out Shelf
Accessible from either side of table; 9" wide.
- B 004 – Drawer
Euro-style drawer with metal sides and laminate bottom and ends. Holds up to 75 lbs.
- *005 – Special Height
36" maximum, 21" minimum. When ordering, state desired height to top of pad on table.
- *006 – Special Length
78" maximum, 48" minimum length. When ordering state desired length (table top).
- *C 008 – Armboard
Adjustable in height and angle. Includes two mounting brackets.
- D 009 – Nosehole
Permits patients to lie face down without turning head. 9" long x 2" wide.
- *E 010 – Stirrups
Adjustable in height, angle, and length. Heels are coated; stirrup bars are chrome plated. Specify type of table when ordering.

- F 020 – Small Pillow
Dimensions: 14" x 12" x 3"
- 060 – Large Pillow
Dimensions: 14" x 21" x 3"
- G 030 – Paper Dispenser
Accommodates paper rolls from 18" to 21" wide.
- H 040 – Paper Cutter
Clear plastic with elastic straps and buckle
- I 077 – Safety Strap
2" x 92" Strap secures with quick fasten hook and loop fastener.
- J 088 – Flat Top Option
Flat table top option can not ship FedEx® Ground and MUST ship by common carrier with all other normal shipping charges
- * K 086 – Red Oak Option (Classic Series Only)
- * L 079 – Triple Bolted Legs Option (Classic Series Only)
- * M 018 – Removable Mat Top Option
2" thick removable mat is held in place with Velcro® strips (Only available with non adjustable backrest tops)

*option not available on Alpha Series tables.

Upholstery Colors	Base Laminate Colors
 3SB - Slate Blue	 1NT - * Natural
 3DT - Desert Tan	 1DC - * Dark Cherry
 3BK - Black	 1MP - Maple
 3PG - Purplegray	
 3VP - Viscaya Palm	
 3AS - Allspice	
 3BK - Black	
 3SP - Sapphire	
 3GP - Graphite	
 3AB - Aubergine	
 3AL - Alabaster	
 3NL - Neutral	

CLINTON VINYL UPHOLSTERY

Clinton uses only premium, woven knit-backed upholstery, that resists stretching and bunching associated with ordinary vinyl upholstery. Plus, Clinton upholstery is topcoated to resist stains, scuffs and abrasions. Standard flammability specifications meet CAL 117. For additional flammability requirements, call for quote.

ORDERING INFORMATION

- PICK A TOP COLOR**
Upholstery for table or laminate for cabinet
- CHOOSE A LAMINATE FRAME COLOR**
Legs, frame and table backs
- SELECT DOOR & DRAWER LAMINATE COLOR**
NOTE: If one color is selected, it will be used on the entire base
- SELECT DOOR & DRAWER PULLS**
NOTE: If no selection is noted on the order, Gray pulls will be used on solid colors, and Silver Wire pulls will be used on wood tones.



Color Coordinating Upholstery	Color Coordinating Laminates
 3GM - Gunmetal	 2SG - Slate Gray
 3DG - Dove Gray	 2GR - Gray
 3CM - Cream	 2FS - Fossil
 3CS - Clamshell	 2SD - Sand
 3WW - Wedgewood	 2AB - Aztec Blue
 3RB - Royal Blue	 2DB - Deep Blue
 3EM - Emerald	 2IG - Imperial Green
 3MB - Mulberry	 2WR - Wild Rose
 3BG - Burgundy	 2SGA - Sangria

Due to the limitations of the printing process, actual colors may vary. Actual upholstery and laminate samples are available.

Quick Ship Log—3 to 5 days
Quick ship items do have similar item quantity restrictions and also have some option limitations. Call for more information.