

Clinton-Eco is one of the best “Green” lines of medical equipment in the industry. Our tables exceeds the strict 2009 California CARB I mandates by using materials containing no added urea-formaldehyde which along with other eco-friendly features improve indoor air quality. Overall, Clinton provides more green features more options and offers one of the broadest selection of eco-friendly equipment in the industry.

Clinton-Eco Wood Tables feature:

- Clinton exceeds 2009 CARB I standard by using 2011 CARB II compliant plywood and laminated panels with no added urea-formaldehyde
- Water-based or UV-cured topcoats
- PVC-free upholstery
- PVC-free edge tape for shelving

Clinton-Eco Metal Tables feature:

- 100% recyclable metal table frame
- Eco-friendly applied VOC-free powder-coated metal
- PVC-free upholstery
- PVC-free edge tape for shelving
- Clinton exceeds 2009 CARB I standard by using 2011 CARB II compliant laminated panels with no added urea-formaldehyde

Clinton-Eco Cabinet Tables feature:

- Clinton exceeds 2009 CARB I standard by using 2011 CARB II compliant plywoods and laminated panels with no added urea-formaldehyde
- PVC-free upholstery
- PVC-free edge tape for shelving
- 100% recyclable metal components

NOTE: Clinton-Eco upholstery and laminates are available for other products. Please call Clinton for quotes.

CLINTON -ECO BASE FINISHES



Natural

Dark Cherry



Maple

Gray (not available for wood tables)

CLINTON -ECO PVC-FREE UPHOLSTERY



Amethyst

South Pacific



Peacock

Sahara



River Rock

Rum



Saffron

Ink



CLINTON™ *ECO* 

81010



81013



81020



81030



Eco-Friendly Wood Treatment Table

	Length	Height	Width
81010-24	72"	31"	24"
81010-27	72"	31"	27"
81010-30	72"	31"	30"

- Features all hardwood legs, H-bracing and adjustable backrest
- Seamless rounded top corners
- 400 lbs. load capacity under normal use
- Also available with non-adjustable backrest (Option 088)*

Eco-Friendly Wood Treatment Table with Drawers

	Length	Height	Width
81013-27	72"	31"	27"
81013-30	72"	31"	30"

- Features all hardwood legs, laminate shelf, two drawers and adjustable backrest
- Seamless rounded top corners
- 400 lbs. load capacity under normal use
- Also available with non-adjustable backrest (Option 088)*

Eco-Friendly Wood Treatment Table with Shelf

	Length	Height	Width
81020-24	72"	31"	24"
81020-27	72"	31"	27"
81020-30	72"	31"	30"

- Features all hardwood legs, laminate shelf and adjustable backrest
- Seamless rounded top corners
- 400 lbs. load capacity under normal use
- Also available with non-adjustable backrest (Option 088)*

Eco-Friendly Wood Treatment Table with Shelving

	Length	Height	Width
81030-24	72"	31"	24"
81030-27	72"	31"	27"
81030-30	72"	31"	30"

- Features all hardwood legs, laminate shelf, storage compartments and adjustable backrest
- Seamless rounded top corners
- 400 lbs. load capacity under normal use
- Also available with non-adjustable backrest (Option 088)*

* 088 Optional non-adjustable top also available. Tables with 088 option MUST ship via common carrier.

CLINTON™ *eco* 

83010



83013



83020



83030



Eco-Friendly Steel Treatment Table

	Length	Height	Width
83010-27	72"	31"	27"
83010-30	72"	31"	30"

- Features powder-coated, all-steel frame, and adjustable backrest
- Seamless rounded top corners
- Paper dispenser included
- 450 lbs. load capacity under normal use
- Also available with non-adjustable backrest (Option 088)*

Eco-Friendly Steel Treatment Table with Drawers

	Length	Height	Width
83013-27	72"	31"	27"
83013-30	72"	31"	30"

- Features powder-coated, all-steel frame, grey laminate shelf, 2 easy-clean drawers and adjustable backrest
- Seamless rounded top corners
- Paper dispenser included
- 450 lbs. load capacity under normal use
- Also available with non-adjustable backrest (Option 088)*

Eco-Friendly Steel Treatment Table with Shelf

	Length	Height	Width
83020-27	72"	31"	27"
83020-30	72"	31"	30"

- Features powder-coated, all-steel frame, gray laminate shelf and adjustable backrest
- Seamless rounded top corners
- Paper dispenser included
- 450 lbs. load capacity under normal use
- Also available with non-adjustable backrest (Option 088)*

Eco-Friendly Steel Treatment Table with Shelving

	Length	Height	Width
83030-27	72"	31"	27"
83030-30	72"	31"	30"

- Features powder-coated, all-steel frame, gray laminate shelf, storage compartments and adjustable backrest
- Seamless rounded top corners
- Paper dispenser included
- 450 lbs. load capacity under normal use
- Also available with non-adjustable backrest (Option 088)*

* 088 Optional non-adjustable top also available.

Tables with 088 option MUST ship via common carrier.

CLINTON™ *ECO* 

89044



Eco-Friendly Cabinet Style Treatment Table

	Length	Height	Width
89044-24	72"	31"	24"
89044-27	72"	31"	27"
89044-30	72"	31"	30"

- Features an adjustable backrest, 2 storage compartments with sliding doors and an adjustable shelf in each compartment
- Seamless rounded top corners
- 425 lbs. load capacity under normal use
- Also available with non-adjustable backrest (Option 088)**

89074



Eco-Friendly Cabinet Style Treatment Table

	Length	Height	Width
89074-24	72"	31"	24"
89074-27	72"	31"	27"
89074-30	72"	31"	30"

- Features 4 doors, adjustable backrest and 2 storage compartments each with an adjustable shelf
- Seamless rounded top corners
- 425 lbs. load capacity under normal use
- Also available with non-adjustable backrest (Option 088)**

89105-30



Eco-Friendly Cabinet Style Treatment Table

	Length	Height	Width
89105-30	72"	31"	30"

- Features 3 drawers, adjustable backrest and 1 storage compartment with 1 adjustable shelf
- Seamless rounded top corners
- 425 lbs. load capacity under normal use
- Also available with non-adjustable backrest (Option 088)**

89110




Eco-Friendly Space Saver Table

89110-24	72" extended/56" folded		
Length	Height	Width	Foot Rest Width
72"	31"	24"	21"
89110-27	72" extended/56" folded		
Length	Height	Width	Foot Rest Width
72"	31"	27"	23"

- Features 1 adjustable shelf, 2 doors, adjustable footrest and adjustable backrest
- Seamless rounded top corners
- Foot rest may be elevated to 3 positions
- 350 lbs. load capacity under normal use
- Also available with non-adjustable backrest (Option 088)**

TABLE OPTIONS



- *A 003 – Slide Out Shelf
Accessible from either side of table; 9" wide.
- B 004 – Drawer
Euro-style drawer with metal sides and laminate bottom and ends. Holds up to 75 lbs.
- *005 – Special Height
36" maximum, 21" minimum. When ordering, state desired height to top of pad on table.
- *006 – Special Length
78" maximum, 48" minimum length. When ordering state desired length (table top).
- *C 008 – Armboard
Adjustable in height and angle. Includes two mounting brackets.
- D 009 – Nosehole
Permits patients to lie face down without turning head. 9" long x 2" wide.
- *E 010 – Stirrups
Adjustable in height, angle, and length. Heels are coated; stirrup bars are chrome plated. Specify type of table when ordering.
- F 020 – Small Pillow
Dimensions: 14" x 12" x 3"
- 060 – Large Pillow
Dimensions: 14" x 21" x 3"
- G 030 – Paper Dispenser
Accommodates paper rolls from 18" to 21" wide.
- H 040 – Paper Cutter
Clear plastic with elastic straps and buckle
- I 077 – Safety Strap
2" x 92" Strap secures with quick fasten hook and loop fastener.
- J 088 – Flat Top Option
Flat table top option can not ship FedEx® Ground and MUST ship by common carrier with all other normal shipping charges
- *K 086 – Red Oak Option (Classic Series Only)
- *L 079 – Triple Bolted Legs Option (Classic Series Only)
- *M 018 – Removable Mat Top Option
2" thick removable mat is held in place with Velcro® strips (Only available with non adjustable backrest tops)

*option not available on Alpha Series tables.