


ETA STYLE LINE SERIES

Contemporary design and easy clean laminate surfaces highlight the ETA (Easy-To-Assemble) Style Line Series of tables. Its flat panel components use bolt-through, pre-threaded hardware for easy and secure assembly. Plus, Clinton's Style Line Series is available in 9 coordinated laminate and vinyl colors to complement any decor.



Pre-mounted heavy-duty, all steel adjustable backrest support system



088 Optional non-adjustable top also available. Tables with 088 option MUST ship via common carrier with normal charges.

9001



Treatment Table with H-Brace

	Length	Height	Width
9001-27	72"	31"	27"
9001-30	72"	31"	30"

- Features all laminate surfaces and adjustable backrest
- Quick Ship in Gray and Natural only
- 350 lbs. load capacity under normal use
- Also available with non-adjustable backrest (Option 088)

9004



Treatment Table with Shelving

	Length	Height	Width
9004-27	72"	31"	27"
9004-30	72"	31"	30"

- Features shelf with storage compartments and adjustable backrest.
- Quick Ship in Gray and Natural only
- 350 lbs. load capacity under normal use
- Also available with non-adjustable backrest (Option 088)

9020



Treatment Table with Shelf

	Length	Height	Width
9020-27	72"	31"	27"
9020-30	72"	31"	30"

- Features full length storage shelf and adjustable backrest
- Quick Ship in Gray and Natural only
- 350 lbs. load capacity under normal use
- Also available with non-adjustable backrest (Option 088)

9013




Treatment Table with Drawers

	Length	Height	Width
9013-27	72"	31"	27"
9013-30	72"	31"	30"

- Features two easy clean, plastic drawers, shelf and adjustable backrest
- Quick Ship in Gray and Natural only
- 350 lbs. load capacity under normal use
- Also available with non-adjustable backrest (Option 088)
- This model ships assembled by common carrier only

TABLE OPTIONS



- *A 003 – Slide Out Shelf
Accessible from either side of table; 9" wide.
- B 004 – Drawer
Euro-style drawer with metal sides and laminate bottom and ends. Holds up to 75 lbs.
- *005 – Special Height
36" maximum, 21" minimum. When ordering, state desired height to top of pad on table.
- *006 – Special Length
78" maximum, 48" minimum length. When ordering state desired length (table top).
- *C 008 – Armboard
Adjustable in height and angle. Includes two mounting brackets.
- D 009 – Nosehole
Permits patients to lie face down without turning head. 9" long x 2" wide.
- *E 010 – Stirrups
Adjustable in height, angle, and length. Heels are coated; stirrup bars are chrome plated. Specify type of table when ordering.
- F 020 – Small Pillow
Dimensions: 14" x 12" x 3"
- 060 – Large Pillow
Dimensions: 14" x 21" x 3"
- G 030 – Paper Dispenser
Accommodates paper rolls from 18" to 21" wide.
- H 040 – Paper Cutter
Clear plastic with elastic straps and buckle
- I 077 – Safety Strap
2" x 92" Strap secures with quick fasten hook and loop fastener.
- J 088 – Flat Top Option
Flat table top option can not ship FedEx® Ground and MUST ship by common carrier with all other normal shipping charges
- *K 086 – Red Oak Option (Classic Series Only)
- *L 079 – Triple Bolted Legs Option (Classic Series Only)
- *M 018 – Removable Mat Top Option
2" thick removable mat is held in place with Velcro® strips (Only available with non adjustable backrest tops)

*option not available on Alpha Series tables.

Upholstery Colors	Base Laminate Colors
 3SB - Slate Blue	 1NT - * Natural
 3DT - Desert Tan	 1DC - * Dark Cherry
 3BK - Black	 1MP - Maple
 3PG - Purplegray	
 3VP - Viscaya Palm	
 3AS - Allspice	
 3BK - Black	
 3SP - Sapphire	
 3GP - Graphite	
 3AB - Aubergine	
 3AL - Alabaster	
 3NL - Neutral	

CLINTON VINYL UPHOLSTERY

Clinton uses only premium, woven knit-backed upholstery, that resists stretching and bunching associated with ordinary vinyl upholstery. Plus, Clinton upholstery is topcoated to resist stains, scuffs and abrasions. Standard flammability specifications meet CAL 117. For additional flammability requirements, call for quote.

ORDERING INFORMATION

- PICK A TOP COLOR**
Upholstery for table or laminate for cabinet
- CHOOSE A LAMINATE FRAME COLOR**
Legs, frame and table backs
- SELECT DOOR & DRAWER LAMINATE COLOR**
NOTE: If one color is selected, it will be used on the entire base
- SELECT DOOR & DRAWER PULLS**
NOTE: If no selection is noted on the order, Gray pulls will be used on solid colors, and Silver Wire pulls will be used on wood tones.



Color Coordinating Upholstery	Color Coordinating Laminates
 3GM - Gunmetal	 2SG - Slate Gray
 3DG - Dove Gray	 2GR - Gray
 3CM - Cream	 2FS - Fossil
 3CS - Clamshell	 2SD - Sand
 3WW - Wedgewood	 2AB - Aztec Blue
 3RB - Royal Blue	 2DB - Deep Blue
 3EM - Emerald	 2IG - Imperial Green
 3MB - Mulberry	 2WR - Wild Rose
 3BG - Burgundy	 2SGA - Sangria

Due to the limitations of the printing process, actual colors may vary. Actual upholstery and laminate samples are available.

Quick Ship Log—3 to 5 days
Quick ship items do have similar item quantity restrictions and also have some option limitations. Call for more information.