

## IMAGING TABLES



### Power Imaging Table with Window Drop

85100	Length	Width	Height Range
Overall	72"	27"	18"-35"

- All-welded steel base, top frame and scissors lift
- Open base design allows easy, close access
- Hands-free foot control
- Built-in hinged window drop. 10" x 10"
- Metal slide bar supports window drop
- Optional locking casters
- 450 lbs. load capacity under normal use

### Power Imaging Table with Window Drop and Adjustable Backrest

85200	Length	Width	Height Range
Overall	72"	27"	18"-35"
Backrest	38"	27"	

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- All welded steel base, top frame and scissors lift
- Hands-free foot control
- Built-in hinged window drop. 10" x 10"
- Metal slide bar supports window drop
- Optional locking casters
- 450 lbs. load capacity under normal use

### Power Imaging Table with 3 Section Top

85300	Length	Width	Height Range
Overall	76"	27"	18"-35"
Drop	13"	27"	
Backrest	38"	27"	

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Drop section, adjustable from either side and raises up to 30° from flat position
- All welded steel base, top frame and scissors lift
- Hands-free foot control
- Built-in hinged window drop. 10" x 10"
- Metal slide bar supports window drop
- Optional locking casters
- 450 lbs. load capacity under normal use

### Multi-Use Power Imaging Table with Stirrups

85309	Length	Width	Height Range
Overall	76"	27"	18"-35"
Drop	13"		
Backrest	38"	27"	

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Drop section, adjustable from either side and raises up to 30° from flat position
- All welded steel base, top frame and scissors lift
- Hands-free foot control
- Built-in hinged window drop. 10" x 10"
- Metal slide bar supports window drop
- Retractable, 2-position, adjustable stirrups
- Optional locking casters
- 450 lbs. load capacity under normal use

## POWER TABLE OPTIONS



080

### Safety Strap

	Length	Width
080	92"	2"

- 2"x 92" strap is secured with quick hook and loop fasteners
- 2 safety straps per set



083

### 2-Hook IV Pole and Mounting Bracket

#### 083

- Chrome pole with height adjustment
- Removable bracket mounts to either side of head end
- Also available with 4-hook top (must specify)

### Hand-Held Height Control

#### 084

- Alternative to foot control
- Stores easily under top when not in use



084



087

### Locking Casters

#### 087

- Heavy duty, rubber wheel, easy-roll casters




082

### Chrome Safety Rail

#### 082

- 3 height position
- Plunger height adjustment

## TABLE OPTIONS



- \*A 003 – Slide Out Shelf  
Accessible from either side of table; 9" wide.
- B 004 – Drawer  
Euro-style drawer with metal sides and laminate bottom and ends. Holds up to 75 lbs.
- \*005 – Special Height  
36" maximum, 21" minimum. When ordering, state desired height to top of pad on table.
- \*006 – Special Length  
78" maximum, 48" minimum length. When ordering state desired length (table top).
- \*C 008 – Armboard  
Adjustable in height and angle. Includes two mounting brackets.
- D 009 – Nosehole  
Permits patients to lie face down without turning head. 9" long x 2" wide.
- \*E 010 – Stirrups  
Adjustable in height, angle, and length. Heels are coated; stirrup bars are chrome plated. Specify type of table when ordering.
- F 020 – Small Pillow  
Dimensions: 14" x 12" x 3"
- 060 – Large Pillow  
Dimensions: 14" x 21" x 3"
- G 030 – Paper Dispenser  
Accommodates paper rolls from 18" to 21" wide.
- H 040 – Paper Cutter  
Clear plastic with elastic straps and buckle
- I 077 – Safety Strap  
2" x 92" Strap secures with quick fasten hook and loop fastener.
- J 088 – Flat Top Option  
Flat table top option can not ship FedEx® Ground and MUST ship by common carrier with all other normal shipping charges
- \*K 086 – Red Oak Option (Classic Series Only)
- \*L 079 – Triple Bolted Legs Option (Classic Series Only)
- \*M 018 – Removable Mat Top Option  
2" thick removable mat is held in place with Velcro® strips (Only available with non adjustable backrest tops)

\*option not available on Alpha Series tables.

Upholstery Colors	Base Laminate Colors
 3SB -Slate Blue	 1NT -* Natural
 3DT-Desert Tan	 1DC-* Dark Cherry
 3BK -Black	 1MP -Maple
 3PG -Purplegray	
 3VP -Viscaya Palm	
 3AS -Allspice	
 3BK -Black	
 3SP -Sapphire	
 3GP -Graphite	
 3AB -Aubergine	
 3AL -Alabaster	
 3NL -Neutral	

## CLINTON VINYL UPHOLSTERY

Clinton uses only premium, woven knit-backed upholstery, that resists stretching and bunching associated with ordinary vinyl upholstery. Plus, Clinton upholstery is topcoated to resist stains, scuffs and abrasions. Standard flammability specifications meet CAL 117. For additional flammability requirements, call for quote.

### ORDERING INFORMATION

- PICK A TOP COLOR**  
Upholstery for table or laminate for cabinet
- CHOOSE A LAMINATE FRAME COLOR**  
Legs, frame and table backs
- SELECT DOOR & DRAWER LAMINATE COLOR**  
NOTE: If one color is selected, it will be used on the entire base
- SELECT DOOR & DRAWER PULLS**  
NOTE: If no selection is noted on the order, Gray pulls will be used on solid colors, and Silver Wire pulls will be used on wood tones.



Color Coordinating Upholstery	Color Coordinating Laminates
 3GM -Gunmetal	 2SG -Slate Gray
 3DG -Dove Gray	 2GR -Gray
 3CM -Cream	 2FS -Fossil
 3CS -Clamshell	 2SD -Sand
 3WW -Wedgewood	 2AB -Aztec Blue
 3RB -Royal Blue	 2DB -Deep Blue
 3EM -Emerald	 2IG -Imperial Green
 3MB -Mulberry	 2WR -Wild Rose
 3BG -Burgundy	 2SGA -Sangria

Due to the limitations of the printing process, actual colors may vary. Actual upholstery and laminate samples are available.

Quick Ship Log—3 to 5 days  
Quick ship items do have similar item quantity restrictions and also have some option limitations. Call for more information.