

OPEN BASE POWER TABLES

FEATURES

Clinton Open Base Power Tables are built tough and engineered for the power you need. With a smooth, dependable scissors lift mechanism, Clinton Open Base Power Tables provide superior performance for a wide variety of medical applications. And they're built here in the USA, so there is never an issue with supply. Plus, Clinton offers them in your choice of 20 upholstery colors to complement any decor or color scheme.

- Ultra-strong uni-frame constructed base and top frame provides a rigid patient platform
- Smooth electric height adjustment moves straight up and down
- Pneumatic, backrest cylinder are internally mounted on the frame for a cleaner appearance
- Hardened bolts and self lubricating pivot points for years of service
- Dependable 110 volt motor
- Optional casters, safety rails and more
- Made in the USA by Clinton Industries
- Low heights help achieve ADA accessibility tax credit
- Dependable scissors lift mechanism that does not need to roll on the floor
- Many different top configurations and options available to fit your needs
- 3-year limited warranty on frame, and 2-years limited warranty on electric components
- Ships fully assembled with no upcharge





Power Table with One Piece Top

80100	Length	Width	Height Range
Overall	72"	27"	18"-35"

- All welded steel base, top frame and scissors lift
- Open base design allows easy close access
- Hands-free foot control
- Optional locking casters
- 450 lbs. load capacity under normal use



Power Table with Adjustable Backrest

80200	Length	Width	Height Range
Overall	72"	27"	18"-35"
Backrest	38"	27"	

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- All welded steel base, top frame and scissors lift
- Hands-free foot control
- Optional locking casters
- 450 lbs. load capacity under normal use



Power Table with Adjustable Backrest & Drop Section

80330	Length	Width	Height Range
Overall	76"	27"	18"-35"
Backrest	38"	27"	
Drop Section	13"	27"	

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Drop section, adjustable from either side and raises up to 30° from flat position
- All welded steel base, top frame and scissors lift
- Nose hole is standard in drop section
- Hands-free foot control
- Optional locking casters
- 450 lbs. load capacity under normal use



Space Saver Power Table with Drop Section

80220	Length	Width	Height Range
Overall	76"	27"	18"-35"
Drop Section	13"	27"	

- Drop section adjustable from either side and raises up to 30° from flat position
- All welded steel base, top frame and scissors lift
- Nose hole is standard in drop section
- Hands-free foot control
- Optional locking casters
- 450 lbs. load capacity under normal use



Power Casting Table with Laminate Leg Rest

80350	Length	Width	Height Range
Overall	76"	27"	18"-35"
Backrest	38"	27"	
Foot Rest	13"	26"	

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Leg section is adjustable from either side and raises up to 30° from flat position
- Easy clean laminate leg rest surface
- Hands-free foot control
- 450 lbs. load capacity under normal use



Power Table with Stirrups, Adjustable Backrest & Drop Section

80360	Length	Width	Height Range
Overall	76"	27"	18"-35"
Backrest	38"	27"	
Foot Rest	13"		

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Drop section is adjustable from either side and raises up to 30° from flat position
- Retractable 2-position adjustable stirrups
- Hands-free foot control
- Optional locking casters
- 450 lbs. load capacity under normal use



Extra Wide Bariatric Straight Top Power Table

84100-34	Length	Width	Height Range
Overall	72"	34"	18"-35"

84100-40	Length	Width	Height Range
Overall	72"	40"	18"-35"

- All welded uni-frame constructed steel base and top frame
- 1" plywood top for extra support
- Hands-free foot control
- Optional locking casters
- 550 lbs. load capacity under normal use



Extra Wide Bariatric Power Table with Adjustable Backrest

84200-34	Length	Width	Height Range
Overall	72"	34"	18"-35"
Backrest	38"	34"	

84200-40	Length	Width	Height Range
Overall	72"	40"	18"-35"
Backrest	38"	40"	

- Dual cylinder, pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- All welded uni-frame constructed steel base and top frame
- 1" plywood top for extra support
- Hands-free foot control
- Optional locking casters
- 550 lbs. load capacity under normal use



Extra Wide Bariatric Power Table with Adjustable Backrest & Drop Section

84430-34	Length	Width	Height Range
Overall	76"	34"	18"-35"
Backrest	38"	34"	
Drop Section	13"	34"	

84430-40	Length	Width	Height Range
Overall	76"	40"	18"-35"
Backrest	38"	40"	
Drop Section	13"	40"	

- Dual cylinder pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Drop section, adjustable from either side and raises up to 30° from flat position
- All welded uni-frame constructed steel base and top frame
- Hands-free foot control
- 1" plywood top for extra support
- Optional locking casters
- 550 lbs. load capacity under normal use

POWER TABLE OPTIONS



Safety Strap

	Length	Width
080	92"	2"

- 2"x 92" strap is secured with quick hook and loop fasteners
- 2 safety straps per set



2-Hook IV Pole and Mounting Bracket

083

- Chrome pole with height adjustment
- Removable bracket mounts to either side of head end
- Also available with 4-hook top (must specify)

Hand-Held Height Control

084

- Alternative to foot control
- Stores easily under top when not in use



Locking Casters

087

- Heavy duty, rubber wheel, easy-roll casters




Chrome Safety Rail

082

- 3 height position
- Plunger height adjustment

TABLE OPTIONS



- *A 003 – Slide Out Shelf
Accessible from either side of table; 9" wide.
- B 004 – Drawer
Euro-style drawer with metal sides and laminate bottom and ends. Holds up to 75 lbs.
- *005 – Special Height
36" maximum, 21" minimum. When ordering, state desired height to top of pad on table.
- *006 – Special Length
78" maximum, 48" minimum length. When ordering state desired length (table top).
- *C 008 – Armboard
Adjustable in height and angle. Includes two mounting brackets.
- D 009 – Nosehole
Permits patients to lie face down without turning head. 9" long x 2" wide.
- *E 010 – Stirrups
Adjustable in height, angle, and length. Heels are coated; stirrup bars are chrome plated. Specify type of table when ordering.
- F 020 – Small Pillow
Dimensions: 14" x 12" x 3"
- 060 – Large Pillow
Dimensions: 14" x 21" x 3"
- G 030 – Paper Dispenser
Accommodates paper rolls from 18" to 21" wide.
- H 040 – Paper Cutter
Clear plastic with elastic straps and buckle
- I 077 – Safety Strap
2" x 92" Strap secures with quick fasten hook and loop fastener.
- J 088 – Flat Top Option
Flat table top option can not ship FedEx® Ground and MUST ship by common carrier with all other normal shipping charges
- *K 086 – Red Oak Option (Classic Series Only)
- *L 079 – Triple Bolted Legs Option (Classic Series Only)
- *M 018 – Removable Mat Top Option
2" thick removable mat is held in place with Velcro® strips (Only available with non adjustable backrest tops)

*option not available on Alpha Series tables.

Upholstery Colors	Base Laminate Colors
 3SB -Slate Blue	 1NT -* Natural
 3DT-Desert Tan	 1DC-* Dark Cherry
 3BK -Black	 1MP -Maple
 3PG -Purplegray	
 3VP -Viscaya Palm	
 3AS -Allspice	
 3SP -Sapphire	
 3GP -Graphite	
 3AB -Aubergine	
 3AL -Alabaster	
 3NL -Neutral	

CLINTON VINYL UPHOLSTERY

Clinton uses only premium, woven knit-backed upholstery, that resists stretching and bunching associated with ordinary vinyl upholstery. Plus, Clinton upholstery is topcoated to resist stains, scuffs and abrasions. Standard flammability specifications meet CAL 117. For additional flammability requirements, call for quote.

ORDERING INFORMATION

- PICK A TOP COLOR**
Upholstery for table or laminate for cabinet
- CHOOSE A LAMINATE FRAME COLOR**
Legs, frame and table backs
- SELECT DOOR & DRAWER LAMINATE COLOR**
NOTE: If one color is selected, it will be used on the entire base
- SELECT DOOR & DRAWER PULLS**
NOTE: If no selection is noted on the order, Gray pulls will be used on solid colors, and Silver Wire pulls will be used on wood tones.



Color Coordinating Upholstery	Color Coordinating Laminates
 3GM -Gunmetal	 2SG -Slate Gray
 3DG -Dove Gray	 2GR -Gray
 3CM -Cream	 2FS -Fossil
 3CS -Clamshell	 2SD -Sand
 3WW -Wedgewood	 2AB -Aztec Blue
 3RB -Royal Blue	 2DB -Deep Blue
 3EM -Emerald	 2IG -Imperial Green
 3MB -Mulberry	 2WR -Wild Rose
 3BG -Burgundy	 2SGA -Sangria

Due to the limitations of the printing process, actual colors may vary. Actual upholstery and laminate samples are available.

Quick Ship Log—3 to 5 days
Quick ship items do have similar item quantity restrictions and also have some option limitations. Call for more information.