

SHROUDED POWER TABLES

Clinton provides plenty of smooth dependable power to our line of Shrouded Power Tables. Engineered and built tough in the USA by Clinton Industries, our tables provide durability, convenience and versatility for a wide variety of applications. And not only do Clinton Shrouded Power Tables perform at superior levels, but they look great too. They are available in a wide variety of upholstery colors to complement any hospital, office or clinic.

FEATURES

- Ultra-strong uni-frame constructed base and top frame provides a rigid patient platform
- Smooth electric height adjustment moves straight up and down
- · Pneumatic, backrest cylinder is internally mounted on the frame for a cleaner appearance
- · Hardened bolts and self lubricating pivot points for years of service
- Dependable 110 volt motor
- · Optional casters, safety rails and more
- Made in the USA by Clinton Industries
- Low heights help achieve ADA accessibility tax credit
- Dependable scissors lift mechanism that does not need to roll on the floor
- Many different top configurations and options available to fit your needs
- 3-year limited warranty on frame and 2-year



* Optional hand-held control available.



DRE, Inc. 1800 Williamson Court Louisville, KY 40223 USA





Shrouded Power Table with One Piece Top

81100	Length	Width	Height Range
Overall	72"	27"	21"-34"

- Textured ABS plastic encloses lift mechanism
- All welded steel base and top frame
- · Hands-free foot control
- Enclosed base makes cleaning easier
- Optional locking casters
- 450 lbs. load capacity under normal use



Shrouded Power Table with Adjustable Backrest

81200	Length	Width	Height Range
Overall	72"	27"	21"-34"
Backrest	38"	27"	

- Textured ABS plastic encloses lift mechanism
- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- All welded steel base and top frame
- Hands-free foot control
- Enclosed base makes cleaning easier
- Optional locking casters
- 450 lbs. load capacity under normal use



Shrouded Power Table with Adjustable Backrest & Drop Section

81330	Length	Width	Height Range
Overall	76"	27"	21"-34"
Backrest	38"	27"	
Drop Section	13″	27"	

- Textured ABS plastic encloses lift mechanism
- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Drop section, adjustable from either side and raises up to 30° from flat position
- All welded steel base and top frame
- Nose hole is standard in drop section
- Hands-free foot control
- Enclosed base makes cleaning easier
- Optional locking casters
- 450 lbs. load capacity under normal use







Shrouded Space Saver Power Table with Drop Section

81220	Length	Width	Height Range
Overall	76"	27"	21"-34"
Drop Section	13"	27"	

- Textured ABS plastic encloses lift mechanism
- Drop section adjustable from either side and raises up to 30° from flat position
- All welded steel base and top frame
- Nose hole is standard in drop section
- Hands-free foot control
- · Enclosed base makes cleaning easier
- Optional locking casters
- 450 lbs. load capacity under normal use



Shrouded Power Adjustable Casting Table with Laminate Surface Leg Rest

81350	Length	Width	Height Range
Overall	76"	27"	21"-34"
Backrest	38"	27"	
Foot Rest	13″	26"	

- Textured ABS plastic encloses lift mechanism
- •Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- All welded steel base and top frame
- •Leg section is adjustable from either side and raises up to 30° from flat position
- Easy clean laminate leg rest surface
- Hands-free foot control
- · Enclosed base makes cleaning easier
- Optional locking casters
- 450 lbs. load capacity under normal use



Shrouded Power Table with Stirrups, Adjustable Backrest & Drop Section

81360	Length	Width	Height Range
Overall	76"	27"	21"-34"
Backrest	38"	27"	
Foot Rest	13"		

- Textured ABS plastic encloses lift mechanism
- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Drop section is adjustable from either side and raises up to 30° from flat position
- Retractable 2-position adjustable stirrups
- Hands-free foot control
- Optional locking casters
- 450 lbs. load capacity under normal use



POWER TABLE OPTIONS













Safety Strap

	Length	Width
080	92″	2"

- 2"x 92" strap is secured with quick hook and loop fasteners
- 2 safety straps per set

2-Hook IV Pole and Mounting Bracket

083

- Chrome pole with height adjustment
- Removable bracket mounts to either side of head end
- Also available with 4-hook top (must specify)

Hand-Held Height Control

084

- Alternative to foot control
- Stores easily under top when not in use

Locking Casters

087

· Heavy duty, rubber wheel, easy-roll casters


Chrome Safety Rail

082

- 3 height position
- Plunger height adjustment







*A 003 - Slide Out Shelf

Accessible from either side of table; 9" wide.

B 004 – Drawer

Euro-style drawer with metal sides and laminate bottom and ends. Holds up to 75 lbs.

*005 – Specialleight

36" maximum, 21" minimum. When ordering, state desired height to top of pad on table.

* 006 – Special Length

78" maximum, 48" minimum length. When ordering state desired length (table top).

*C 008 - Armboard

Adjustable in height and angle. Includes two mounting brackets.

D 009 – Nosehole

Permits patients to lie face down without turning head. 9" long x 2" wide.

*E 010 - Stirrups

Adjustable in height, angle, and length. Heels are coated; stirrup bars are chrome plated. Specify type of table when ordering.

F 020 - Small Pillow

Dimensions: 14" x 12" x 3"

060 – Large Pillow

Dimensions: 14" x 21" x 3"

G 030 – Paper Dispenser

Accommodates paper rolls from 18" to 21" wide.

H 040 – Paper Cutter

Clear plastic with elastic straps and buckle

I 077 – Safety Strap

2" x 92" Strap secures with quick fasten hook and loop fastener.

J 088 – Flat Top Option

Flat table top option can not ship FedEx* Ground and MUST ship by common carrier with all other normal shipping charges

- * K 086 Red Oak Option (Classic Series Only)
- * L 079 Triple Bolted Legs Option (Classic Series Only)

* M 018 - Removable Mat Top Option

2" thick removable mat is held in place with Velcro* strips (Only available with non adjustable backrest tops)

*option not available on Alpha Series tables.



Color CLINTON VINYL UPHOLSTERY Upholstery **Base Laminate** Coordinating Coordinating Colors Colors Upholstery Laminates Clinton uses only premium, woven knit-backed upholstery, that resists stretching and bunching associated with ordinary vinyl upholstery. Plus, Clinton upholstery is topcoated to 3SB -Slate Blue 1NT -* Natural 2SG-Slate Gray 3GM-Gunmetal resist stains, scuffs and abrasions. Standard flammability specifications meet CAL 117. For additional flammability requirements, call for quote. 3DT-Desert Tan 1DC-* Dark Cherry 3DG-Dove Gray 2GR -G ray ORDERING INFORMATION 1. PICK A TOP COLOR Upholstery for table or laminate for cabinet 2. CHOOSE A LAMINATE FRAME COLOR 3BK -Black 1MP -Maple Legs, frame and table backs 3CM -Cream 2FS -Fossil 3. SELECT DOOR & DRAWER LAMINATE COLOR NOTE: If one color is selected, it will be used on the entire base 3PG-Purplegray 4. SELECT DOOR & DRAWER PULLS NOTE: If no selection is noted on the order, 3CS -Clamshell 2S D-Sand Gray pulls will be used on solid colors, and Silver Wire pulls will be used on wood tones. 3VP -Viscaya Palm 3WW -Wedgewood 2AB -Aztec Blue 3AS - Allspice 3RB-Royal Blue 2DB-Deep Blue 3SP -Sapphire 3EM -Emerald 2IG-Imperial Green 3GP-Graphite 3MB -Mulberry 2W R-Wild Rose 3AB - Aubergine 3BG-Burgundy 2SGA-Sangria Due to the limitations of the printing process, 3AL - Alabaster actual colors may vary. Actual upholstery and laminate samples are available. Quick Ship Log—3 to 5 days Quick ship items do have similar item quantity restrictions and also have some option limitations. Call for more information.



3NL -Neutral