

Helps you meet the demands of today's healthcare environment



The Dinamap Compact TS is capable of monitoring vital signs, including non-invasive blood pressure (NIBP), pulse oximetry (SpO₂) and predictive body temperature.

Suitable for adult, pediatric, and neonatal use, the Compact TS offers unparalleled accuracy and sensitivity for monitoring hypertensive patients. This model handles patient movement very well.



Allows you to take all connected vitals with the press of one button

STAT mode for rapid BP readings in emergency situations

Single determinations or automatic monitoring of blood pressure

User-adjustable target inflation pressure

Oral/rectal predictive temperature

Built-in printer available

Powerful, lightweight and easy to use

STAT NIBP allows interval blood pressure monitoring

Display shows last three NIBP determination

Trends store up to 100 past vital sign readings